

ON-GOING PROJECTS R&D AND TECHNICAL SUPPORT

« BITS'N PIECES »

S. Happel

TrisKem Users Group Meeting

14th of September 2010

Chester (UK)

Outline


Results of some on-going projects

- Ga-68 purification
- Zr separation
- Sn separation

Ga purification

- Context: purification of Ga-68 obtained from Ge-68/Ga-68 generator
- Too high levels of Fe and Zn in Ga-68 fraction
- Literature:
 - Cation exchange
 - Anion exchange
 - UTEVA (McAlister/Horwitz)


Ga purification



Elution study 2 mL TEVA column; 6M HCl, Elution condition: 0.1M HCl

- **Good separation from Zn (stays on column)**
- **No separation from Ge, little separation from Fe**
 - Fe separation improvement would lead to Ga loss


Ga purification



Elution study 2 mL UTEVA column; 4M HCl, Elution condition: 0.1M HCl

- Good separation from Ge
- No separation from Zn and Fe


Ga purification



Elution study 2 mL LN column; 5M HCl, Elution condition: 0.5M HCl

- Good separation from Ge, Fe (remains on resin) and Zn
- > 80% Ga in 1 mL 0.5M HCl
- Quantitative elution in 2 mL 0.5M HCl

Ga purification



Elution study 2 mL LN column; 5M HCl, Elution condition: 1M HCl

- Good separation from Ge, Fe and Zn
- > 95% Ga in 1 mL 1M HCl
- Quantitative elution in 2 mL 1M HCl


Zr-93

- Long-lived radionuclide
- Decommissioning and radioactive waste
- LSC or ICP-MS (→ Nb-93!)
- Combination of literature methods
 - Vajda et.al. (UTEVA)
 - Le Fèvre/Pin (UTEVA)
 - Bombard et. al. (TEVA)
 - NAS (AIX)

Zr separation via UTEVA/TEVA

- **1st separation on UTEVA**
 - Advantage: no selectivity for Fe
- **Method following Vajda et. al.**
 - Modified elution
 - Important for loading solution: quantitative F⁻ removal / masking!
- **Load: 5mL 9M HNO₃ / 100 μL H₂O₂**
 - Multi-element solutions approx. 100μg per element
- **Rinse: - 5 mL 9M HNO₃ / 50 μL H₂O₂**
 - 2 x 5 mL 9M HNO₃
 - 5 mL 5M HNO₃
- **Elution: 4 x 5 mL 0.01M HNO₃ / 0.1M oxalic acid**
- **All fractions analysed by ICP-MS**


Zr separation on UTEVA



➤ **Zr elution works fine**


- Near quantitative recovery in 5 mL 0.01M HNO₃ / 0.1M oxalic acid

Zr separation on UTEVA



- No Zr / U separation
- Ag separation not quantitative

Zr separation on UTEVA




- No Zr / Hf separation
 - Not conform with literature....
- Sb(III) separation to be confirmed
- Good Nb separation

Zr separation on TEVA – UTEVA eluate


- **2nd separation on TEVA**
 - Advantage: different selectivity than UTEVA
- **Method following Bombard et. al.**
 - Modified elution
- **Column preconditioning: 8M LiNO₃ / 0.01M HNO₃**
- **Load: 3mL 9M HCl**
 - UTEVA eluate from first test evaporated and converted to 9M HCl
- **Rinse: 2 x 1 mL 9M HCl**
- **Elution: 3 x 5 mL 0.05M HCl**
- **All fractions analysed by ICP-MS**
- **Additionally applied to multi-element solution**

Zr separation on TEVA



- Good separation from Ag
- No separation from U
 - U/Zr separation already on UTEVA
 - Modification of Zr elution from TEVA

Zr separation on TEVA




- Separation of Nb, Mo, Sb and Hf from Zr problematic under chosen conditions
 - Hf and Sb separation should be improved by increase of 9M HCl rinsing step volume
 - Mo behavior not conform with literature – control of oxidation state

Zr separation on UTEVA – modified elution

- No Zr/U separation with 0.01M HNO₃ / 0.1M oxalic acid
- Horwitz et al. show that Zr can be eluted with 2M HNO₃ / 0.02 M oxalic acid while U stays on UTEVA
- Load: 5mL 9M HNO₃ / 100 μL H₂O₂
 - Multi-element solutions approx 100μg per element
- Rinse: - 5 mL 9M HNO₃ / 50 μL H₂O₂
 - 2 x 5 mL 9M HNO₃
 - 5 mL 5M HNO₃
- Elution: 4 x 5 mL 2M HNO₃ / 0.03M oxalic acid
- All fractions analysed by ICP-MS

Zr separation via UTEVA/TEVA




- Good U/Zr separation
- Still no Zr/Hf separation
- Th separation delicate
- Zr elution in 5 mL 2M HNO₃ / 0.03M oxalic acid

Zr separation on TEVA – modified elution

- **NAS describes anion exchange method based on**
 - Loading from 9M HCl
 - Zr Elution in 7M HCl
- **Column preconditioning: 8M LiNO₃ / 0.01M HNO₃**
- **Load: 3mL 9M HCl**
 - Multi-element solutions approx 100µg per element
- **Rinse: 2 x 1 mL 9M HCl**
- **Elution: 4 x 5 mL 7M HCl**
- **All fractions analysed by ICP-MS**

Zr separation on TEVA



- Th and U separation very good
- Interesting Zr/Hf elution behaviour
- Zr recovered in 10 mL 7M HCl
- To be confirmed...

Conclusion

- **Combining UTEVA and TEVA seems promising**
- **Method still under work, needs to be improved**
- **Influence of matrix and HF to be evaluated**
- **Separation on UTEVA:**
 - **Improvement of Zr/Hf separation**
 - **Confirmation of U and Th separation from Zr**
 - **Confirmation of Nb/Zr separation**
 - **Influence of HF and matrix...**
- **Separation on TEVA**
 - **Preconditioning with 9M HCl**
 - **Increased volume of 9M HCl rinsing step**
 - **Confirmation for Mo and Hf elution behaviour**
 - **Zr Elution via 0.05M HCl or 7M HCl**

Sn separation



- Context: Sn-121m in decommissioning samples
- High levels of matrix elements (incl. Ca and Fe)
- Literature:
 - Ikaza/Yata: TBP resin
 - Dissertation M. Langer: TRU resin
 - Load: 5mL 8M HCl (Fe and Sn retained, Co, Cd and Ni pass)
 - Rinse: 25 mL 8M HCl (removes remaining Co, Cd and Ni)
 - Rinse: 2.5M HCl (Fe elution)
 - Sn Elution with 20 mL 0.1M HCl

Sn separation – 1st Elution study

- Preparation of TBP resin
- 2 mL TBP resin column
- Preconditioned with 5 mL 2M HCl
- 5 mL sample solution containing (V, Cr, Mn, Fe, Ni, Co, Cu, Zn, Y, Ag, Cd, Sn, Sb, Cs, Ba, Nd, Sm, Eu, Pb and U - each approx 50µg) in 2M HCl
- Sample loaded onto TBP resin column – eluate collected
- Rinsing with 15 mL 2M HCl (collected in 1 mL steps)
- Elution with 20 mL 0.1M HCl (collected in 1 mL steps)



- All fractions measured by ICP-MS

Results – 1st Elution study



➤ Overall good selectivity for Sn, elution in 5 mL 0.1M HCl possible

Results – 1st Elution study



- No Fe/Sn separation
- Zn separation to be improved....

Sn - Fe / Zn separation

- Ghersini/Braun: In 0.1M H₂SO₄ TBP has very low selectivity for Fe and Zn, and high selectivity for Sn
- Potential method:
 - 2 mL column TBP resin
 - Preconditioned with 5 mL 2M HCl
 - 5 mL sample loading solution 2M HCl
 - Rinsing with 2 x 5 mL 2M HCl
 - Rinsing with 5 mL 0.1M H₂SO₄
 - Elution with 5 mL 0.1M HCl
- Test currently under way
- For iron/metal samples: load from 0.1M H₂SO₄, followed by rinse with 2M HCl to avoid saturation of column with Fe?
- Further testing of method
- Testing of TRU resin as described by M Langer

Very special thanks!!

- R. Streng
- C. Dirks
- A. Zulauf

at University Marburg