

Rapid Method for Po-210 in Urban Matrices

Sherrod L. Maxwell, PhD

Co-author Dan McAlister, PhD

Triskem User's Group Meeting 2018

September 21, 2018

Recent News- Potential RDD Plans

- “Brussels attacks: Belgium fears Isis seeking to make ‘dirty’ nuclear bomb” March-2016

<http://www.independent.co.uk/news/world/europe/brussels-attacks-belgium-fears-isis-seeking-to-make-dirty-nuclear-bomb-a6951661.html>

- “An official at Belgium’s Federal Agency for Nuclear Control told *The Times*: **“When you start filming someone in the way they did, the logical conclusion is that they wanted to abduct that person and to obtain radioactive material.”**



Need for Rapid Methods

- 1998: Chechen authorities foiled a possible terrorist act when they found and defused a mine attached to a container "full of radioactive substances" near Grozny,
 - terrorists stole radioactive waste from Russian RADON storage facility
- ***"Mexico Finds Stolen Radioactive Material Amid Dirty Bomb Fears"***
<http://www.voanews.com/content/nuclear-material-stolen-in-mexico/1803195.html>, 12/4/13
"A dirty bomb detonated in a major city could cause mass panic, as well as serious economic and environmental consequences"
- ***Are we ready??***



List of Threat Radionuclides

Alpha Emitters			Beta/Gamma Emitters		
Radionuclide	Half-Life	Emission Type	Radionuclide	Half-Life	Emission Type
²⁴¹ Am	432.6 y	α, γ	Ac-227 [†]	21.77 y	β, γ
²⁵² Cf	2.64 y	α, γ			
²⁴² Cm	163 d	α	Ce-141 [*]	32.51 d	β, γ
²⁴³ Cm	29.1 y	α, γ	Ce-144 [‡]	284.9 d	β, γ
²⁴⁴ Cm	18.10 y	α	Co-57 [*]	271.7 d	ε, γ, x-ray
²³⁷ Np	2.14×10 ⁶ y	α, γ, x-ray	Co-60 [*]	5.271 y	β, γ
²¹⁰ Po [*]	138.4 d	α	Cs-134 [*]	2.065 y	β, γ
²³⁸ Pu	87.7 y	α	Cs-137 [§]	30.07 y	β, γ
²³⁹ Pu	2.41×10 ⁴ y	α	H-3 [*]	12.32 y	β only
²⁴⁰ Pu	6.56×10 ³ y	α	I-125 [*]	59.40 d	ε, γ, x-ray
²²⁶ Ra [†]	1.60×10 ³ y	α, γ	I-129 [†]	1.57×10 ⁷ y	β, γ, x-ray
²²⁸ Th	1.912 y	α, γ	I-131 [*]	8.021 d	β, γ
²³⁰ Th	7.538×10 ⁴ y	α, γ	Ir-192 [*]	73.83 d	β, γ
²³² Th	1.405×10 ¹⁰ y	α	Mo-99 [†]	65.94 h	β, γ
²³⁴ U	2.455×10 ⁵ y	α	P-32 [*]	14.26 d	β only
²³⁵ U	7.038×10 ⁸ y	α, γ	Pd-103 [*]	16.99 d	β, γ
²³⁸ U	4.468×10 ⁹ y	α	Pu-241	14.29 y	α, β
U-Nat	—	α	Ra-228 [†]	5.75 y	β only
			Ru-103 [†]	39.26 d	β, γ
			Ru-106 [†]	373.6 d	β only
			Se-75 [*]	119.8 d	ε, γ
			Sr-89 [*]	50.53 d	β only
			Sr-90 [†]	28.79 y	β only
			Tc-99 [*]	2.11×10 ⁵ y	β only

The half-lives of the nuclides are given in years (y), days (d) or hours (h)

* No radioactive progeny or progeny not analytically useful.

† Radioactive progeny with short half-lives, and the progeny may be used as part of the detection method for the parent.

* Radioactive progeny not used for quantification, only screening.

§ Radioactive progeny used for quantification only, not screening.

Po-210

- ^{210}Po ($T_{1/2}=138.38$ days)
 - naturally-occurring alpha emitting radionuclide used in industrial applications
 - **static elimination**
 - ^{210}Po is present in the environment as a result of the ^{238}U decay chain
 - one of the most toxic naturally-occurring radionuclides
 - an *ingestion* intake as small as $1\ \mu\text{g}$ of *Po-210* may be *lethal* for the most radiosensitive members of the population



Risk of Po-210 Attack

- **Former Russian intelligence officer Alexander Litvinenko died in a London hospital on November 23, 2006**
 - *Deliberately poisoned with ^{210}Po*
 - *Police discovered those involved in this crime had spread ^{210}Po over many locations in London*
 - *Incident caused widespread psychological concern and a heightened political and public health response*
- **Naturally-occurring radionuclides such as ^{226}Ra and ^{210}Po can also be used in a radiological dispersive device (RDD) as part of a terrorist attack**



From a Tea Pot....

- “Pinpoint Homicide”
- ~1000 people requested follow-up (urine analysis)
- 137 people contaminated (above background)
- 19 Environmental Teams, 24/7 for weeks
- 47 sites evaluated for contamination
- Many locations “sealed” for ~ 5 years
- Some people view this event as a “Micro-Terrorist” event
- **Can we imagine the impact of a Po-210 RDD?**
 - A non-explosive dispersive attack is also a serious concern


Current Analytical Method Challenges

- What are the challenges with Po-210 assay?
 - Acid digestion of solid samples- Po volatilization
- Methods typically involve spontaneous auto-deposition of ^{210}Po onto silver or other metal disks followed by alpha spectrometry
 - auto-deposition times range from 90 minutes to 24 hours or more, at times with yields that may be less than desirable
 - If sample interferences are present, decreased yields and degraded alpha spectra can occur, due to an unpredictable increase in thickness in the deposited layer
- Problematic, especially during a radiological emergency,
 - where rapid detection of ^{210}Po is critical and the safety of the public is at stake
 - method rigor and minimizing rework are extremely important

- International Atomic Energy Agency (IAEA) ALMERA (Analytical Laboratories for the Measurement of Environmental Radioactivity) network
 - *administered proficiency testing for the rapid analysis of ^{210}Po in water after the Litvinenko event.*
 - *ALMERA network, established by the IAEA in 1995, is a technical collaboration of existing institutions and makes available to Member States*
 - *a worldwide network of analytical laboratories capable of providing reliable and timely analysis of environmental samples in the event of an accidental or intentional release of radioactivity.*

ALMERA Network PT Study: How did they do?

Labs in PT study



- ***Only 70 % of the results reported by the 33 labs met the IAEA requirements***
- ***Only 19 of the 33 labs met the one week target***
- *22 of 33 participants reported results which “fit the purpose of rapid detection of Po-210 in water”*
- ***Reporting times were usually about 10 days for most labs, but some labs took over 20 days and still did not report acceptable results***

Challenge: Polonium Volatilization

“Polonium is often lost in the initial digestion because it has a relatively low volatilization temperature, 180 C”

Martin A, Blanchard RL (1969) The thermal volatilisation of caesium-137, Po-210 and lead-210 from in vivo labelled samples. *Analyst* 94(1119):441–446. doi:10.1039/AN9699400441

Mabuchi H (1963) The volatility of some Po compounds. *J Inorg Nucl Chem* 25:657–660. doi:10.1016/0022-1902(63)80154-2

What can we do??

Savannah River Environmental Laboratory Approach

- **Alkaline fusion**
 - Fast, rugged, **inhibits Po volatilization**
 - Quick preconcentration –Fe(OH)₃
- **DGA Resin- k' in HCl for Po ~100,000**
- **Po-210 on Sr Resin**
 - eluted with 6M HNO₃, autodeposition, 50-70% yields,
 - **k' in HCl using DGA Resin much higher**
 - Vajda, N., La Rosa, J., Zeisler, R., Danesi, P., Kis-Benedek, G.Y., (1997) A Novel technique for the simultaneous determination of ²¹⁰Pb and ²¹⁰Po using a crown ether, J. Environ. Radioactive. 37: 355-372
- ***Use fast microprecipitation instead of autodeposition***
 - Can be eluted from DGA Resin and collected quickly (BiPO₄)
 - Optimized in this work

Use of DGA Resin and BiPO_4 Microprecipitation

J Radioanal Nucl Chem
DOI 10.1007/s10967-013-2644-2

Rapid determination of ^{210}Po in water samples

Sherrod L. Maxwell • Brian K. Culligan •
Jay B. Hutchison • Robin C. Utsey •
Daniel R. McAlister

Calcium phosphate ppt. and DGA Resin
+ BiPO_4 microprecipitation

Received: 22 May 2013
© Akadémiai Kiadó, Budapest, Hungary 2013

Challenges

- **SRS solid sample fusion methods**
 - Preconcentration of actinides with Fe/Ti hydroxide
 - Fe/Ti/silicates removal with LaF_3
 - Fe can interfere with some resin methods
 - *Fe (III) has synergistic effect that may cause Ca(II) ions to be retained more than normal*
 - Colloidal Si can be a nightmare
- **Cannot precipitate Po as a fluoride to remove Fe/silicates!**
- **How do we overcome this problem??**

Potential Solution


- **Fluoride can stabilize silicates in solution**
- **Can we add fluoride to the sample load solution**
 - to keep Si in solution and eliminate Si gel/colloids and clogging
 - *and still retain Po on DGA?*
- **Test a dilute HCl-HF solution to load Po^{+4} onto DGA Resin**
 - *Po (IV) retention is very high even in low HCl*
 - *If Po complexing is stronger with chloride...may be ok*

Success!

- *F⁻ does not interfere significantly with Po⁺⁴ retention on DGA*
- **Silicates pass through column without colloidal solids/gel formation!**
- *Fe has minimal retention on DGA Resin in low HCl*
 - *Fe(OH)₃ can be employed*
 - *Does not interfere*




Fe on DGA Resin in HCl



D.R. McAlister and E.P. Horwitz, "Synergistic enhancement of the extraction of trivalent lanthanides and actinides by tetra(n-octyl)diglycolamide from chloride media," *Solv. Extr. Ion Exch.*, 26(1), 12-24 (2008).

New Po-210 Fusion Method




Maintain alkalinity
on heat ramp to avoid Po loss

Keep HCl low

Use F⁻ to keep Si 'happy'
in solution

Po Retention on DGA Resin





We can elute Po
with 0.05M HNO_3
however 10-15% would
stay on resin

How to optimize?


*Add trace oxalate to
0.05M HNO_3*

Bi Retention on DGA Resin



Fast Po-210 Purification using DGA Resin

1M HCl-Boric Acid rinse
to eliminate any traces
of ThF₄ on Resin



Now using HNO₃-oxalate to enhance Po elution from DGA

Added indicator to ensure optimal pH


....Bi elutes as well

Po Elution- DGA Resin

Elution of Po-210 on 2 mL TODGA Cartridge				
		Po-210 LSC		
Fraction	total mL	%		
Load*	25	0.47		
Rinse 1M HCl - 0.25M Boric Acid	5	0.05		
Rinse 0.25M HCl	20	0.14		
Rinse 0.5M HCl - 0.01M HF	10	0.14		
Rinse 8M HNO ₃	10	0.02		
Strip 0.05M HNO ₃ - 0.005M Biocalate	5	96.0	total 20 mL strip 96.7	
Strip 0.05M HNO ₃ - 0.005M Biocalate	5	0.3		
Strip 0.05M HNO ₃ - 0.005M Biocalate	5	0.2		
Strip 0.05M HNO ₃ - 0.005M Biocalate	5	0.2		
cartridge	2	2.4		
Precondition cartridge 10 mL 1M HCl				
*Load = 5 mL 4M HCl-0.2M HF + 20 mL 0.5M HCl - 0.02M HF				
DGA lot DNSR16A				

85-94% recovery with 0.05M HNO₃ only


Po Elution Curve on DGA Resin



A = Precondition 10 mL 1M HCl. B = Load 5 mL 4M HCl – 0.02M HF + 23 mL 0.5M HCl – 0.02M HF. C = 5 mL 1M HCl. D = 20 mL 0.25M HCl. E = 5 mL 0.5M HCl – 0.02M HF. F = 10 mL 8M HNO₃. G = 20 mL 0.05M HNO₃ – 0.005M ammonium bioxalate.

Will Bi Isotopes Interfere?

- Bi-210 ($T_{1/2} = 5.1$ days) 100% beta emission abundance
- Bi-214 ($T_{1/2} = 19.9$ min) 0.02% alpha emission abundance
- Po-214 ($T_{1/2} = 136.6$ μ sec) 99.99% alpha abundance, 7.69 MeV
- Bi-212 ($T_{1/2} = 60.55$ min) 35.94% alpha abundance, ~6.05 MeV
- Po-212 ($T_{1/2} = 2.98$ E-7s) 100% alpha abundance, ~8.78 MeV



Results for Po-210 Spiked in Concrete Samples

Sample ID	²⁰⁹ Po Yield (%)	²¹⁰ Po Reference Value (mBq g ⁻¹)	²¹⁰ Po Measured Value (mBq g ⁻¹)	²¹⁰ Po Native Value (mBq g ⁻¹)	²¹⁰ Po Corrected Value (mBq g ⁻¹)	Difference (%)
1	99.6	366.7	486.9	110.9	376.0	2.5
2	99.8	366.7	491.0	110.9	380.1	3.7
3	86.4	366.7	485.1	110.9	374.2	2.0
4	91.8	366.7	494.7	110.9	383.8	4.7
5	81.3	366.7	495.8	110.9	384.9	5.0
6	78.9	366.7	476.2	110.9	365.3	-0.4
7	88.8	366.7	529.1	110.9	418.2	14.0
8	79.8	366.7	495.8	110.9	384.9	5.0
9	94.0	366.7	484.7	110.9	373.8	1.9
10	88.8	366.7	481.0	110.9	370.1	0.9
11	83.6	366.7	462.1	110.9	351.2	-4.2
12	96.6	366.7	465.1	110.9	354.2	-3.4
13	90.4	366.7	465.1	110.9	354.2	-3.4
14	82.8	366.7	439.9	110.9	329.0	-10.3
15	86.5	366.7	460.7	110.9	349.8	-4.6
Avg	88.6		480.9			0.89
SD	6.8		20.9			5.69
% RSD	7.8		4.3			

Results for Po-210 Spiked in Granite Samples

Sample	²⁰⁹ Po Yield	²¹⁰ Po Reference Value	²¹⁰ Po Measured Value	²¹⁰ Po Native Value	²¹⁰ Po Corrected Value	Difference
ID	(%)	(mBq g ⁻¹)	(mBq g ⁻¹)	(mBq g ⁻¹)	(mBq g ⁻¹)	(%)
1	96.4	366.7	460.7	97.49	363.2	-1.0
2	92.2	366.7	478.0	97.49	380.5	3.8
3	104.7	366.7	464.4	97.49	366.9	0.1
4	92.8	366.7	494.3	97.49	396.8	8.2
5	98.8	366.7	480.0	97.49	382.5	4.3
6	88.5	366.7	462.1	97.49	364.6	-0.6
Avg	95.6		473.3			2.47
SD	5.7		13.2			3.60

Results for Po-210 Spiked in Soil Samples


Sample	²⁰⁹ Po Yield	²¹⁰ Po Reference Value	²¹⁰ Po Measured Value	²¹⁰ Po Native Value	²¹⁰ Po Corrected Value	Difference
ID	(%)	(mBq g ⁻¹)	(mBq g ⁻¹)	(mBq g ⁻¹)	(mBq g ⁻¹)	(%)
1	84.2	366.7	417.7	64.57	353.1	-3.7
2	99.3	366.7	432.2	64.57	367.6	0.3
3	74.0	366.7	442.9	64.57	378.3	3.2
4	78.0	366.7	433.6	64.57	369.0	0.6
5	78.2	366.7	461.0	64.57	396.4	8.1
Avg	82.7		437.5			1.69
SD	10.0		15.9			4.35

Po-210 Spiked In Soil-Enhanced Precipitation

Sample ID	²¹⁰ Po Yield (%)	²¹⁰ Po Reference Value (mBq g ⁻¹)	²¹⁰ Po Measured Value (mBq g ⁻¹)	²¹⁰ Po Native Value (mBq g ⁻¹)	²¹⁰ Po Corrected Value (mBq g ⁻¹)	Difference (%)
1	88.5	366.7	441.0	64.57	376.4	2.7
2	98.7	366.7	410.7	58.46	352.2	-3.9
3	90.0	366.7	431.8	58.46	373.3	1.8
4	94.0	366.7	414.4	58.46	355.9	-2.9
Avg	92.8		424.5			-0.60
SD	4.6		14.4			3.32

Added 25 mg Ca at Fe(OH)₃ precipitation step....

Po-210 Alpha Spectra Granite Samples



Analytical Method Timeline

- | | |
|--|------------------|
| • Spiking and Alkaline Fusion | 45 minutes |
| • Preconcentration as $\text{Fe}(\text{OH})_3$ | 45 minutes |
| • DGA Resin Separation | 90 minutes |
| • BiPO_4 microprecipitation | 30 minutes |
| • Sample preparation time | 3 hr. 30 minutes |
| • Alpha Spectrometry counting | 60-1000 minutes |



Summary

- **New fast method for Po-210 in urban matrices**
 - Rugged alkaline fusion to minimize Po volatilization
 - *Concrete*
 - *Granite*
 - *Soil*
 - Fast alpha source preparation
 - *Eliminates potential autodeposition issues, can save time, improve yields, spectra*
 - *BiPO₄ microprecipitation*
 - *Potential to stack Sr Resin and DGA Resin for Pb-210*
- **Future work**
 - Limestone, Marble
 - Glass fiber air filters