

# Characterization of a TBP Resin and development of methods for the separation of actinides and the purification of Sn

Dirks C.<sup>a</sup>, Vajda N.<sup>b</sup>, E. Kovács-Széles<sup>c</sup>, A. Bombard<sup>d</sup>, S. Happel<sup>d</sup>

<sup>a</sup> Philipps-University, Marburg, GERMANY, <sup>b</sup> Radanal Kft, Budapest, HUNGARY, <sup>c</sup> Centre for Energy Research, Hungarian Academy of Sciences, Budapest, HUNGARY, <sup>d</sup> TrisKem International, Bruz, FRANCE

Philipps




Universität  
Marburg



- General
- Resin characterisation
- Application I – Pu in drinking water
- Application II – Sn separation

# General

- TBP used in PUREX process (LLX)
- TRU Resin (contains TBP) used for Sn separation in geology/archeology
  - Elimination of matrix elements and isobaric interferences
- Determination of long-lived Sn isotopes in rad waste
  - Focus on matrix removal and elimination of isobaric interferences
- Sn-117m separation for use in nuclear medicine
  - Focus on Sn/Cd/Sb separation





TBP (Tributyl Phosphate)

# Resin characterisation

- Determination of  $D_w$  values for various elements
  - Multi-element solutions (HCl and HNO<sub>3</sub>) for ICP-MS
    - 10 µg/mL of each element: Al, As, B, Ba, Ca, Cd, Co, Cr<sup>3+</sup>, Cs, Cu, Fe, Ga, Li, Mg, Mn, Na, Ni, Pb, Rb, Sr, Th, U, V, Zr
    - 10 µg/mL of each element: B, Ge, Mo, Nb, P, Re, S, Si, Ta, Ti, W, Zr
  - Pu(IV), Np(IV), Th(IV) and U(VI) via LSC
  - 50mg resin contacted with 1.3 or 1.5mL solution for ≥ 1h
  - Centrifugation and filtration
  - Dilution with H<sub>2</sub>O for ICP-MS measurements (multi-element solutions),
  - Evaporation and dilution in 0.1M HNO<sub>3</sub> for LSC measurements (Pu, Np, Th and U)

# Resin characterisation


➤  $D_w$  values of the actinides in  $\text{HNO}_3$  and  $\text{HCl}$



- $D_w$  Pu > 100 in  $\geq 6\text{M HNO}_3/\text{HCl}$
- $D_w$  U < 100 all acids/acid concentrations

# Resin characterisation

## ➤ Anionic interferences



## ➤ U in 8M $HNO_3$ :


➤ No/little interference from oxalate, interference from sulfates and especially phosphates

## ➤ Pu in 9M HCl/0,01M $NaNO_2$ :

➤ Interferences, but  $D_w(Pu) > 500 \Rightarrow$  little impact on Pu retention


# Resin characterisation

- $\text{HNO}_3$ : only elements with  $D_w > 10$  shown
- $D_w(\text{Ag}) \sim 500$  in 0.1M  $\text{HNO}_3$  and  $D_w(\text{Pu}) \sim 100$  in 8M  $\text{HNO}_3$
- Other elements show very little affinity in  $\text{HNO}_3$



# Resin characterisation

## ➤ $D_w$ values in various other conditions



# Resin characterisation

## ➤ U capacity

- Batch experiments
- $U_{\text{nat}}$ , 50mg resin, 8M  $\text{HNO}_3$
- Filtration, evaporation, dissolution in 0.1M  $\text{HNO}_3$
- Mixed with LSC cocktail (ProSafe+)
- LSC measurement

## ➤ Results:


- capacity = 71-76 mg U/g dry resin in 8M  $\text{HNO}_3$ 
  - comparable (although significantly lower) to UTEVA resin

# Elution study


- 2mL columns, gravity flow
- Pu-239 (Pu(IV)), U-233 and Th-230
- Presence of Fe(III)
- All eluates collected, evaporated and redissolved in 0.1M HNO<sub>3</sub>
- Mixed with LSC cocktail (ProSafe+)
- Elution procedure:
  - Preconditioning: 10 mL 8M HNO<sub>3</sub>/0.01M NaNO<sub>2</sub>,
  - Loading solution (~ 10 mL)
  - Rinse: 2 x 10 mL 8M HNO<sub>3</sub>/0.01M NaNO<sub>2</sub>,
  - Elution Th: 10 mL 9M HCl/0.01M NaNO<sub>2</sub>,
  - Elution Th+U: 30 mL 9M HCl/0.01M NaNO<sub>2</sub>,
  - Elution Pu: 20 mL 1M HCl

# Elution study

## Separation of U and Th on TBP resin



## Separation of Pu from actinides on TBP resin



# Application I

## ➤ Determination of Pu in drinking water

### ➤ Pre-treatment:

- 300-500mL water acidified with 2.5mL 1M HNO<sub>3</sub> and spiked with Pu-239, Am-241, Th-230 and U-233 (each 2Bq)
- 0,5g Mohr's salt
- 100μL N<sub>2</sub>H<sub>4</sub> (Fe reduction)
- Heating under stirring for 1h
- Fe(OH)<sub>2</sub> precipitation with 2.5mL NH<sub>3</sub> 25% (pH 7)
- Heating, settling over night
- Filtration 0.45μm
- Redissolve in 8M HCl
- Evaporation, redissolve in 3x1mL conc. HNO<sub>3</sub>
- 350μL N<sub>2</sub>H<sub>4</sub> + 1 drop conc. HNO<sub>3</sub> (Test SCN<sup>-</sup>: Fe<sup>2+</sup>)

# Application I


- Determination of Pu in drinking water
  - Pre-treatment:
    - Add 1.5mL conc.  $\text{HNO}_3$  to adjust to 4M  $\text{HNO}_3$  and heat to boiling ( $\text{N}_2\text{H}_4$  is destroyed) and
    - Add 1 drop of conc.  $\text{HNO}_3$  and allow to cool down to RT
    - Addition of  $\text{NaNO}_2$  and conc.  $\text{HNO}_3$  to adjust to 10mL 8M  $\text{HNO}_3$ /0,1M  $\text{NaNO}_2$ .
  - Separation as described
- Results :
  - Chemical yield for Pu ~ 69%.
  - U contamination in Pu source is <1.4%.
  - No Am or Th found
  - Procedure can be performed in 1 day.

# Application I

- Limits of the procedure:
  - With respect to  $D_W(\text{Pu})$  in 8M  $\text{HNO}_3$  loading volume has to be <20 mL.
  - Very iron rich (Fe content > 2 g) and/or large samples (>1-2 g of soil) can probably not be treated on standard 2 mL columns.
  - Pu(IV) valence adjustment time consuming

# Application II

- Sn separation:
  - Method development based on  $D_W$  values obtained via batch experiments
- Results:



**L1:** 4mL loading solution:  
2M HCl/0,01M formic acid

**R1:** 5mL 2M HCl/0,01M formic acid

**R2:** 2mL 2M HCl/0,01M formic acid

**S1:** 3mL 1M HCl


**E1:** 3mL 1M HCl

**E2:** 3mL 1M HCl

**E3:** 4mL 0,1M HCl


**E4:** 2mL 0,1M HCl

# Proposed Sn separation procedure



# Application II

- Main isobaric interference for Sn-126 (ICP-MS determination): Te-126
- Decon factor study to verify Sn/Te separation



➤ Te decon.factor in Sn fraction > 1000

# Application II


- Results

- TBP Resin can be used for the purification of Sn
- Most elements are eluted during load and first rinse (Cd, As, Ag, Ge, Zn, In, ~70% Sb) => 11mL
- Fe/Ga are removed with 9mL 1M HCl
- For Fe rich sample loading under reducing conditions might be necessary
- >90% Sn eluted in 6mL 0,1M HCl
- ~ 30% Sb co-eluted with Sn => control of Sb oxidation state
- Clean Se/Te separation
- On-going project on Sn-126 determination in rad waste via ICP-MS
  - First step AIX, followed by TBP

# Conclusions

- TBP resin characterized with respect to  $D_W$  values of various elements in  $\text{HNO}_3$  and  $\text{HCl}$ , and the influence of selected interferents
- TBP resin is well suited for Sn separation. Optimization of Sn/Sb separation through control of Sb oxidation state (important for work with Sb targets).
- Possibility to separate and determine Pu in drinking water

Thank you for  
your attention!



**TRISKEM**  
SHARING INNOVATION

TRISKEM INTERNATIONAL

Parc de Lormandière Bât. C - Rue Maryse Bastié - Campus de Ker Lann - 35170 Bruz - FRANCE

Tel +33 (0)2.99.05.00.09 - Fax +33 (0)2.99.05.07.27 - [www.triskem-international.com](http://www.triskem-international.com) - email : [contact@triskem.fr](mailto:contact@triskem.fr)