





# Outline

---

- Why Copper?
- Batch Experiments
  - Selectivity
  - Interferences with Ni or Zn
  - Kinetics
- Column Experiments
  - Simulated Ni target
  - Irradiated Ni target
  - Column size
  - Decontamination factors
  - Conversion of Cu resin eluate via AIX
- Summary
- Outlook

# Copper Radionuclides

Isotope	Half Life	Radiation	Source
$^{60}\text{Cu}$	20 min	$\beta^+$ (93%), EC (7%)	cyclotron
$^{61}\text{Cu}$	3.3 hours	$\beta^+$ (62%), EC (38%)	cyclotron
$^{62}\text{Cu}$	9.74 min	$\beta^+$ (98%), EC (2%)	generator/cyclotron
<b><math>^{64}\text{Cu}</math></b>	<b>12.7 hours</b>	<b><math>\beta^+</math>(18%), EC (41%), <math>\beta^-</math> (37%)</b>	<b>reactor/cyclotron</b>
$^{66}\text{Cu}$	5.2 min	$\beta^-$ (100%)	reactor/cyclotron
$^{67}\text{Cu}$	62 hours	$\beta^-$ (100%)	reactor/cyclotron

# Radiopharmaceutical application<sup>[1,2]</sup>

---

- Suitable half-life (  $t_{1/2} = 12.7 \text{ h}$  )
- Multiple decay mode
- Well established coordination chemistry
- **PET imaging**
- **targeted therapy (→ Cu-67)**

# Comparison of production routes of Cu-64<sup>[3,4]</sup>

Nuclear process	Optimum energy range (MeV)	Thick target yield (MBq/μAh)
<b><math>^{64}\text{Ni}(p,n)^{64}\text{Cu}</math></b> <sup>a</sup>	<b>12 → 8</b>	<b>304</b>
$^{64}\text{Ni}(d,2n)^{64}\text{Cu}$ <sup>a</sup>	17 → 11	430
$^{68}\text{Zn}(p,\alpha n)^{64}\text{Cu}$ <sup>a</sup>	30 → 21 <sup>b</sup>	116
$^{66}\text{Zn}(p,2pn)^{64}\text{Cu}$ <sup>a</sup>	52 → 37	316
$^{64}\text{Zn}(d,2p)^{64}\text{Cu}$ <sup>a</sup>	20 → 10	27.1
$^{66}\text{Zn}(d,\alpha)^{64}\text{Cu}$ <sup>a</sup>	13 → 5	13.8
$^{\text{nat}}\text{Zn}(d,x)^{64}\text{Cu}$	25 → 10 <sup>c</sup>	57.0

a: Using highly enriched target material, low enrichment will lead to impurities;

b: Below threshold of  $^{67}\text{Cu}$  impurity *via* the  $^{68}\text{Zn}(p,2p)^{67}\text{Cu}$  reaction;

c: Below thresholds of  $^{61}\text{Cu}$  and  $^{67}\text{Cu}$  impurities *via* the  $^{64}\text{Zn}(d,\alpha n)^{61}\text{Cu}$  and  $^{68}\text{Zn}(d,2pn)^{67}\text{Cu}$  reactions, respectively.

# General procedure for batch experiments

---

## Weight distribution ratios $D_w$

- Weigh approx. 50 mg of the new Cu resin in an 2 ml Eppendorf tube
  - Add 400  $\mu$ l of the acid (e.g. HCl pH 2)
  - Close cap and shake for 30 minutes
  - Add 1ml of the standard solution (e.g. 1 mL Multi-element solution)
  - Close cap and shake for another 30 minutes
  - Withdraw 1 ml of the supernatant, analyze (ICP-MS)
- All distribution factors were determined in triplicate

# Weight distribution ratio

---

$$D_W = \frac{N_{A0} - N_A}{N_A} \times \frac{V}{m_R}$$

- **high  $D_w$  = Extraction**
- **low  $D_w$  = Elution**


$N_{A0}$  = net count rate in the  $A_0$  sample

$N_A$  = net count rate in the sample

$V$  = Volume of the aqueous phase (1.4 mL)

$m_R$  = amount of the resin in g

# $D_w$ coefficients for multi 2 in HCl (each element with 10 $\mu$ g/ml)



## Loading:

- Good selectivity for Cu
- No selectivity for other selected elements

## Elution:

- Low  $D_w$  at low pH
- Tests with 4, 6 and 8 M HCl

Figure 1:  $D_w$  of Cu and selected elements on Cu resin in HCl in varying pH values



# $D_w$ coefficients for multi 3 in $H_2SO_4$ (each element with $10\mu g/ml$ )





Figure 2:  $D_w$  of Cu and selected elements on Cu resin in  $H_2SO_4$  in varying pH values

# $D_w$ coefficients for multi 1 in $\text{HNO}_3$ (each element with $10\mu\text{g}/\text{ml}$ )



- Cu: pH >2 high  $D_w$
- overall good Cu capacity

Figure 1:  $D_w$  of Cu and selected elements on Cu resin in  $\text{HNO}_3$  in varying pH values



# Conclusions I

---

- Stable, high Cu  $D_w$  values at pH > 2
- High selectivity for Cu
- No selectivity for Zn or Ni
- Loading solution: HCl pH 2 (or higher)
- Elution with 6 or 8 M HCl

# $D_w$ Cu – interferences - HCl pH 2





Figure 2:  $D_w$  of Cu on Cu resin in HCl, pH 2 in presence of various amounts of Ni


Figure 3:  $D_w$  of Cu on Cu resin in HCl, pH 2 in presence of various amounts of Zn

# Interferences of Zn or Ni

---

- Even for high amounts of Ni or Zn  $D_w(\text{Cu}) > 1000$  at HCl pH 2
- Up to 1 g target material pro g Resin
  - **negligible interference!**

# Kinetics of the Cu-resin



**fast kinetic**

- 5 min 89 %
- 15 min 97 %
- 60 min 99.8 %

Figure 4: kinetics of Cu on Cu resin in HCl, pH 2

# Elution method

**Loading : 5 mL simulated or dissolved irradiated target solution**

**Rinsing : 2 x 10 mL HCl pH 2**

**Elution Cu : 2x 5 mL HCl 6 or 8 M**




Figure 5: Scheme of elution conditions of the elution study

# Elution study - simulated Ni target




Figure 6: Elution profile obtained for simulated Ni target, 2 mL Cu resin column, elution conditions as described in figure 5




# Elution study - irradiated Ni target

---

- Irradiation of a Ni foil (10 mg) :
  - ( $\varnothing = 13$  mm, 0.025 mm thickness)
  - Cyclotron BC1710 at Jülich
  - $E_p = 15$  MeV
  - 1 h; 0.5  $\mu$ A
- Additionally 170 mg nonirradiated Ni-foil

# Elution study – irradiated Ni target



Almost 100 % Cu in E1

## Loading:

L1 : 5 mL HCl pH 2

## Rinsing:

R1 : 10 mL HCl pH 2


R2 : 10 mL HCl pH 2

## Elution:


E1 : 5 mL HCl 8 M

E2 : 5 mL HCl 8 M


Figure 7 : Elution study: % yield of the elution study with an irradiated Ni target



**Eluate from loading : L 1  
5 mL pH 2 HCl**



**Rinsing : R1**  
**10 mL HCl pH 2**



**Elution : E1**  
**5 mL HCl 8M**

# Conclusions II

---

- ✓ High Cu selectivity
- ✓ Negligible interference of Ni or Zn
- ✓ Fast kinetics
- ✓ Quantative recovery of Cu in elution studies
- ✓ **Very good Cu separation**
- ✓ Ni could be recovered in a small volume of HCl pH 2
- Cu elution volume still rather large

# Optimization of Cu elution volume





Figure 8/9 : Elution study: elution yields in %, simulated Ni target, varying Cu resin amounts

- Use of smaller columns allows reducing Cu elution volume
- Quantitative Cu recovery in 1 mL 8M HCl
- No impact on purity

# Optimization of Cu elution volume – vacuum-assisted flow





Figure 10/11 : Elution study: elution yields in %, simulated Ni and Zn target, 300 mg Cu resin

- Flow rate: 1 mL/min
- 0.3g columns
- Simulated Ni and Zn targets
- > 90% recovery in 1 mL 8M HCl, quantitative recovery in 1.5 mL
- No impact on purity



# Decontamination factors $D_f$

---

- Flow rate: 1 mL/min
- 0.3g columns
- Loading solution: elevated amounts of Ni, Zn, Co, Ga and Au in 5 mL HCl pH 2
- Separation as described before
- ICP-MS measurement
- Calculation of decontamination factors  $D_f$  for Cu fractions
  - Fraction E1 (0.5 mL 8M HCl):
    - $D_f$ : Ni, Co & Zn > 20 000
    - $D_f$ : Au & Ga > 10 000
  - Fraction E2 (0.5 mL 8M HCl):
    - $D_f$ : Ni > 20 000, Co > 40 000, Zn > 70 000, Au > 50 000, Ga > 10 000

# Conversion of Cu eluate

---

- Aim: recovery of Cu in small volume of dilute HCl, water or NaCl solution
- Anion exchange resins (AIX) shows necessary selectivity
- Cu eluate (1mL 8M HCl) from 300 mg Cu resin column directly loaded onto small AIX column
- Rinse with 1 mL 8M HCl
- Elution with 2 x 1 mL H<sub>2</sub>O

# Conversion of Cu eluate via AIX





Figure 12/13 : Elution study: elution yields in %, simulated Ni target, 300 mg Cu resin, varying amounts of AIX resin

- Quantitative load of Cu on AIX from 8M HCl , Cu elution in 1 mL H<sub>2</sub>O
  - Conversion works well
- 150 mg AIX column shows Cu breakthrough in rinsing step
  - Use of 200 mg AIX columns or more advisable
- Next steps:
  - Smaller elution values possible?
  - Additional purification by AIX step?

# Elution study - AIX



- 400 mg AIX
- Cu elution in 0.6 – 0.8 mL water
- Anion exchange conversion step gives additional decontamination from Ni, Zn, Au<sup>18</sup>

# Optimized method

---

- Vacuum-assisted flow (1 – 3 mL/min)
- Cu resin (300 mg)
  - Load from 1 - 2 mL HCl pH2
  - Rinse with 5 mL and 3 mL HCl pH 2
  - Load and rinse contain ~100% Ni (important for Ni-64 recovery)
  - Cu elution in 1 – 1.5 mL 8M HCl
- Conversion on AIX
  - Load from 1 – 1.5 mL 8M HCl
  - Rinse with 2 x 0.5 mL 8M HCl
  - Cu elution in 0.6 – 0.8 mL water (or saline solution)
- Cu yield > 95%, high decontamination factors
- Separation time: 18 minutes

# Short view to the future

---

- Cu-DOTA labeling
- Organic impurities content of the Cu fraction
- Recovery of Ni for target preparation
- Optimisation of flowrates
- Ni-64 and/or Zn target irradiation
- Analytical application (concentration and purification of Cu for ICP-MS)

# Literature

---

- [1] Qaim, S. M.: Decay data and production yields of some nonstandard positron emitters used in PET. *Q. J. Nucl. Med. Mol. Imaging* 52, 111–120 (2008).
- [2] Smith, Suzanne: Molecular imaging with copper-64. *J. of inorganic Biochemistry* 98, 1874–1901 (2004)
- [3] Aslam, M.N., Sudár, S., Hussain, M., Malik, A.A., Shah, H.A. and Qaim, S.M.: Charged particle induced reaction cross section data for production of the emerging medically important positron emitter  $^{64}\text{Cu}$ : A comprehensive evaluation. *Radiochim. Acta* 97, 669–686 (2009)
- [4] Szelecsényi, F., Blessing, G., Qaim, S.M.: Excitation functions of proton induced nuclear reactions on enriched  $^{61}\text{Ni}$  and  $^{64}\text{Ni}$ : Possibility of production of no-carrier-added  $^{61}\text{Cu}$  and  $^{64}\text{Cu}$  at a small cyclotron. *Appl. Radiat. Isot.* 44, 575–580 (1993).
- [5] Horwitz, E.P., McAlister, D.R., Bond, A.H., Barrans, R.E. and Williamson, J.M.: A process for the separation of  $^{177}\text{Lu}$  from neutron irradiated  $^{176}\text{Yb}$  targets. *Appl. Radiat. Isot.* 63(1), 23-36 (2005)
- [6] **Dirks, C., Scholten, B., Happel, S., Zulauf, A., Bombard, A., Jungclas, H.:**  
**Characterisation of a Cu selective resin and its application to the production of Cu-64**  
**J. Radioanal Nucl. Chem. in press**