

Determination of ^{36}Cl in decommissioning samples using a Pyrolyser furnace and extraction chromatographic separations

P E Warwick¹, A Zulauf², S Happel³, I W Croudace¹


¹ GAU Radioanalytical, Southampton, UK

² Philipps Universitat, Marburg, Germany

³ Triskem Ltd, France

Origin of ^{36}Cl

- ^{36}Cl is predominantly produced via neutron activation of naturally occurring ^{35}Cl .
- ^{36}Cl is a long lived (3.02×10^5 y) beta emitting radionuclide ($E_{\text{max}} = 709.6$ keV).
- ^{36}Cl is present in nuclear graphite, concretes, ion exchange resins & desiccants.
- Characterisation of ^{36}Cl in nuclear wastes is important due to its mobility in the geosphere and high soil – plant transfer factor




Current approaches

- Alkaline digestion followed by chemical separation.
- Acid digestion and volatilisation of ^{36}Cl .
- Thermal decomposition of the sample (with or without modifiers) and liberation of Cl species as HCl or Cl_2 .
- Final measurement of ^{36}Cl by liquid scintillation counting (high counting efficiency).
- Potential interferences from ^3H , ^{14}C , ^{35}S , ^{129}I .

Proposed separation

- Thermal decomposition of the sample and desorption of Cl species.
- Samples heated to 900°C in a tube furnace using a controlled heating cycle (10°C/min to 900°C and held for 60 minutes – total time 2 hrs).
- Moist air used to flush system. Sample moistened with 1ml water.
- Combustion products trapped in an alkaline trap. Bubbler connected to furnace tube using ground glass joints.
- ^{36}Cl isolated from other radionuclides using Triskem 'Cl resin' conditioned with AgNO_3
- Total volume of eluent mixed with scintillation cocktail for liquid scintillation analysis.


Thermal desorption of ^{36}Cl



Chlorine-36 standard
as NaCl

Bubbler (20ml 6mM
 Na_2CO_3) changed at
100°C intervals

Recovery of standards



Results for replicate standard analyses


Mean recovery = 86%
(all in 1st bubbler)

Blanks run between samples

Carry over between samples < 0.1%

Residual activity in sample boat ~ 0.3%

Trapping of Cl species



- ³⁶Cl liberated as either Cl₂ or HCl.
- Trapped in 20ml 6 mM Na₂CO₃ (bubblers 1+2) or 1M NaOH (bubbler 3)
- Air flow rate 200 ml/min

Characterisation of the Cl resin

Retention of ^{36}Cl and ^{129}I in 1M H_2SO_4

Isotope	D_w retention
Cl-36	1600
I-129	1980

- quantitative uptake of both isotopes by silver loaded Cl-resin

D_w values for different KSCN concentrations

	Cl-36	I-129
KSCN conc.	D_w elution	D_w elution
0.01M	1.7	12000
0.05M	0.4	15000
0.1M	0.7	4000
0.2M	0.4	9000

- ^{36}Cl is eluted quantitatively at any KSCN concentration
- ^{129}I remains on the resin at any KSCN concentration

D_w values for different Na_2S concentrations

Na_2S conc	Mean D_w
0.04M	40
0.09M	15
0.18M	0.7
0.35M	0.8



- ^{129}I is eluted at elevated Na_2S concentrations

Chloride loading capacity




Analyte	Theoretical value	Experimental value
I-	14.9mg	16.3±1.6mg
Cl-	4.2mg	4.3 ±0.2mg

Loadings are dependent on the quantity of Ag initially loaded onto the resin.
Above values are based on 13mg Ag loading

Separation of ^{36}Cl and ^{129}I from 1M NaOH




Separation from 6 mM Na_2CO_3




Cl recovery 98%

Decontamination factors (I)



Decontamination factors (II)




Additional 0.1M H_2SO_4 wash stage added to remove ^{14}C

Decontamination factors

	³⁶Cl fraction	¹²⁹I fraction
³HTO	> 500	> 2000
¹⁴CO₃	7	5000
¹⁴C modified wash	700	
³⁵S modified wash	1500	1000
³⁶Cl		> 2000
¹²⁹I	1300	

Liquid scintillation analysis

- 5ml 0.1M NH_4SCN eluent mixed with 15ml Proflow P cocktail.
- Some cocktails contain additives which will reduce any residual Ag compounds co-eluted to elemental Ag resulting in a black solution unsuitable for liquid scintillation counting.
- Proflow P (Meridian) mixes well with the NH_4SCN eluent without reducing residual Ag^+ .



Limit of detection


Sample mass (g)	1.0 g
Recovery %	86 %
Background CPM	11 cpm
Efficiency %	98 %
Count time (mins)	180 mins
LOD (Currie)	0.02 Bq/g

Spiked ion exchange resin


Sample type	Expected value	Measured value
Ion exchange resin	4.1 kBq	4.3 ± 0.1 kBq

Analysis of decommissioning samples

- Method tested using desiccant from a reactor site.



LSC spectra for bubbler solution
(dominantly ^3H)



LSC spectra for purified
bubbler solution
(^{36}Cl)

Benefits over existing approaches

- Rapid procedure with analysis requiring 1 day for 6 samples.
- The total bubbler volume is used in the final measurement increasing sensitivity.
- Applicable to decommissioning samples containing other volatile radionuclides.
- Potential for sequential separation and quantification of other volatile radionuclides.

Summary

- Initial studies indicate that ^{36}Cl is effectively liberated from solid matrices using thermal desorption.
- Other volatile radionuclides co-trapped with ^{36}Cl can be efficiently separated using Cl-resin. The resin will also isolate ^{129}I .
- The combination of combustion and isolation of Cl using Cl-columns provides a rapid approach for the separation and purification of ^{36}Cl from solid matrices avoiding the need for time-consuming digestion procedures.
- The Cl and I fractions arising from the separation are readily miscible with commercially available liquid scintillation cocktails.
- Further studies are required to validate the technique and confirm that ^{36}Cl can be quantitatively extracted from the range of materials routinely analysed for ^{36}Cl .

