

Rapid methods for the determination of actinides, radiostrontium and radium in environmental and bioassay samples

Scope

- Actinides and Sr in soil, food, fecal, concrete and brick samples
- Determination of radiostrontium in sea water samples
- Determination of Ra-226 in environmental samples

Actinides and Sr in soil, food, concrete and brick samples

SL Maxwell, BK Culligan, A Kelsey-Wall, PJ Shaw: Rapid radiochemical method for determination of actinides in emergency concrete and brick samples. *Anal Chim Acta.*, 701(1):2011;112-118.


SL Maxwell, BK Culligan, A Kelsey-Wall, PJ Shaw: Rapid determination of actinides in emergency food samples, *J. Radioanal. Nucl. Chem.*, 292(1), 2011, 339-347

S. L. Maxwell and B.K. Culligan: Rapid Method for Determination of Actinides in Fecal Samples, 31/10/12, 58th Annual RRMCC, Fort Collins, CO October 29 to November 2, 2012

- Methods can be adjusted for larger sample masses
- Addition of internal standards and Sr carrier (or Sr-85)
- Mineralisation in furnace at 700°C
- NaOH fusion
- Two co-precipitations for matrix removal
 - $\text{Fe}(\text{OH})_3$ / Ca-Phosphate
 - LaF_3 under reducing conditions ($\text{TiCl}_3 \rightarrow \text{U}(\text{IV})$)
- Dissolution in 3M HNO_3 / 1M $\text{Al}(\text{NO}_3)_3$ / 0.25M boric acid
- Redox (Pu(IV)): $\text{Fe}(\text{II})$ / NaNO_2

- Vacuumbox system
- Stacked TEVA, TRU and DGA Kartuschen -> actinide retention
- Rinse with 3M HNO₃
- Separation of the cartridges (TEVA and TRU/DGA)
 - Th, Pu (Np) purification via TEVA
 - Am/U purification via TRU/DGA
- Microprecipitation
- Eluates from sample load and first rinse (all cartridges) united and evaporated to dryness
- Sr purification on 3 mL Sr Resin column or cartridges (2 mL + 1 mL)

Sample preparation



S. Maxwell et al. 2011

< 1 day

Rapid Fecal Sample Furnace Heating



Place fecal sample in 1L beaker lined with vellum paper

Add Tracers (2x since sample is split later)
[²⁴²Pu/²³⁶Pu, ²⁴³Am, ²³²U]

Place in furnace and heat @ 250°C for 30 min

Ramp 1 to 350°C and heat for 20min *
Ramp 2 to 450°C and heat for 20min
Ramp 3 to 550°C and heat for 45min

* Routine analysis heating option to 550°C in 1 liter overnight without ceramic crucible



Wet ash residue in 1L beaker

Transfer ash to 250mL ceramic crucible



Wet ash residue with 15.8M HNO₃ and 30wt% H₂O₂ on hot plate

Place ceramic crucible in furnace @ 850°C for 1-1.5hr until ash is white

Transfer wet ashed residue (after evaporation to small volume) into 250mL Zr crucible with 15.8M HNO₃

Transfer ash from ceramic crucible to 250mL Zr crucible with 15.8M HNO₃




250mL Zr crucible
Evaporate to dryness

Place Zr crucible in furnace at ~450°C and increase heat (if necessary) until solids are white/light colored



Separation scheme (Sr optional)



Method performance (MAPEP 18 samples)

- Good agreement (bias $15\% \leq B \leq -15\%$)
- High yields for actinides, good yields for Sr

Sample Code	Am yield (%)	Pu yield (%)	U yield (%)	Sr yield (%)
MAPEP-18 soil	96.2±6.33	102.2±10.5	84.0±5.64	60.0±2.8
MAPEP-20	na	na	na	66.0 +/- 6.0
10g baby food	84.6±7.5	93.5±8.1	77.9±13.1	na
10g apple	93.4±9.1	97.5±12.1	88.9±10.9	na
10g squash	88.5±3.5	97.5±5.9	77.9±13.1	na
MAPEP-18 spiked concrete	85.3±6.5	89.6±7.9	76.9±4.4	na
MAPEP-18 spiked brick	93.7±2.9	94.7±9.0	88.1±5.4	na
NRIP fecal	82.7±3.9	96.4±8.2	62.5±7.2	na

S. Maxwell, 2010/11

- Results in < 1d – 2d
- Method can be adapted to ICP-MS


Sr in sea water samples

Maxwell S L, Culligan B K, Utsey R C: Rapid method for the determination of Radiostrontium in sea water samples, 31/10/12, 58th Annual RRMCM, Fort Collins, CO
October 29 to November 2, 2012

Radiostrontium in sea water



- Sea water: 7 – 8 mg Sr / L, 400 mg Ca / L
- ICP-MS for yield
- Preferably samples > 1L for low detection limits
- Preconcentration by coprecipitation
- 2 options for separation:
 - Sr-89/90: combined Sr/DGA resins
 - Sr-90 only: DGA resin
- Measurement via GPC (LSC or Cerenkov also possible)

Rapid Sr-89/90 Sample Preparation Method for Seawater



Maxwell et al. 2012


Rapid Column Separation Method



Load at 1 drop/sec on vacuum box

Shows stacked option

If present, ⁹¹Y is a problem!
If agree, confirmatory but better to allow ingrowth, then separate Y-90




Results Sr-89/90 option

- 1L spiked sea water sample (7,66 mg Sr.L⁻¹), 4 mL Sr resin:
 - 2h counting time
 - Measurement via Sr-90: Yield: 88,8% (+/- 5,9%, N = 11), Bias: 1,2%
 - Measurement via Y-90: Yield: 95,0% (+/- 1,6%, N = 11), Bias: 3,1%
 - Good correspondance
- 2L sea water sample (7,70 mg Sr.L⁻¹), 6 mL Sr resin:
 - 2h counting time
 - Yield: 81,9% (+/- 5,0%, N = 4), Bias: 4,2%
- Measurement via GPC
- MDAs:
 - 1L sea water (2 x 2 mL cartridges)
 - MDAs: 9.1 mBq.L⁻¹ (2h count), 4,4 mBq.L⁻¹ (8h count), 3.0 mB.L⁻¹ (1000 min count)
 - 2L sea water (3 x 2 mL cartridges)
 - MDAs: 9.1 mBq.L⁻¹ (2h count), 4,4 mBq.L⁻¹ (8h count), 3.0 mB.L⁻¹ (1000 min count)
 - 6L sea water (three 2L aliquots combined after purification)
 - MDAs: 1.5 mBq.L⁻¹ (8h count), 1.0 mB.L⁻¹ (1000 min count)

Sr-89/90 option

- Similar methods suggested for environmental water samples
 - Groska J, Molnar Z, Bokori E, Vajda N: Simultaneous determination of ^{89}Sr and ^{90}Sr : comparison of methods and calculation techniques, Journal of Radioanalytical and Nuclear Chemistry, March 2012, Volume 291 (3),707-715
 - T. O'Brien et al.: The rapid determination of Strontium-89 and Strontium-90 in Environmental Samples. Presented at the MARC IX conferences, Kailua-Kona, USA, 29/03/12
- Measurement by Cerenkov counting possible
 - Sr-89 and Y-90 via Cerenkov
 - Very low interference of Sr-90 on Sr-89
 - Advantageous in case of high Sr-89/90 activity ratios

Rapid Column Separation for ^{90}Sr (^{90}Y) – DGA Only Option



Sr-90 (Y-90) DGA resin only option

- **1 to 10 liter method (DGA Resin only)**

- 2 liter aliquot requires one 2 ml DGA Resin cartridge

- *MDA with GFPC and 120 minute count = 9.1 mBq/L*
- *MDA with GFPC and 480 minute count = 4.4 mBq/L*
- *MDA with GFPC and 1000 minute count = 3.0 mBq/L*

- 10 liter aliquot (5 x 2 liter aliquots combined after purification)

- *MDA with GFPC and 480 minute count = 0.88 mBq/L*
- *MDA with GFPC and 1000 minute count = 0.61 mBq/L*

- *< 1mBq/L ⁹⁰Sr with 10L seawater aliquot and < 6 hour sample preparation*


Maxwell et al. 2012


Rapid determination of Ra-226 in environmental samples

- For solid samples use of MnO_2 resin not possible
 - High matrix load after sample dissolution, precipitation at pH 7
- Solid samples frequently contain elevated amounts of Ba
 - Problematic for preparation of source for α -spectrometry
 - Polyatomic Interferences at ICP-MS measurements
- Ba removal necessary
 - Ba/Ra separation (e.g. SR Resin)
 - Ba-133 can not be used as internal standard
 - Alternative: Ra-225/At-217 (from Th-229), advantage: α -Spectrometry

Rapid determination of Ra-226 in environmental samples

- Rapid method Sherrod Maxwell (SRS)
 - Filter, 5g vegetation, 1g soil, brick or concrete, 150 mL water samples
 - Ashing (2h 700°C, wet ashing, 5 – 10 min 600°C)
 - NaOH fusion in Zr crucible
 - Carbonate precipitation
 - Cation exchange (Ca removal)
 - Optional: SR Resin (for Ba-rich samples)
 - LN Resin (Ac, Ca,... removal)
 - Microprecipitation and α -Spectrometry






Maxwell, 2012

Results spiked real samples

Matrix	Chemical yield / %	Obtained result / mBq per sample	Reference value / mBq per sample	Bias to ref. value / %
Vegetables	87.1 (5.7)	72.8 (5.1)	73.8	-1.2
Concrete	84.6 (6.8)	180.6 (8.0)	184.5	-2.1
Brick	86.5 (6.6)	77.8 (4.6)	73.8	5.5
Air filter	76.7 (4.2)	77.1 (6.2)	73.8	4.5
Soil	75.3 (1.9)	184.9 (6.2)	184.5	0.2
Water	91.8 (6.7)	70.9 (3.7)	73.8	-3.9

Maxwell, 2012

- Yields between 75 and 90%
- Good agreement with reference values
- Clean spectra



Maxwell, 2012

Спасибо за внимание!

Вопросы?



www.triskem-international.com



http://www.linkedin.com/company/triskem-international?trk=hb_tab_compy_id_2897456

TRISKEM INTERNATIONAL

Parc de Lormandière - Bât.C - Rue Maryse Bastié - Campus de Ker Lann - 35170 BRUZ - FRANCE

Tél. +33(0)2 99 05 00 09 - Fax. +33(0)2 99 05 07 27 - www.triskem-international.fr - email : contact@triskem.fr