

Update on resins and methods for the purification and QC of radionuclides for use in diagnostics and therapy

2025 Australasian Cyclotron Users Group
Heidelberg (Australia)

Steffen Happel
22/05/2025

**2025 Australasian
Cyclotron Users Group**



Triskem International



- Based in Rennes (France)
- Independent company since 2007 (before part of Eichrom)
- Main product line: extraction chromatographic resins.
 - Also producing other separation materials
- Staff : 25
- R&D and TechSupport group:
 - 3 RadChem PhD, 2 Technicians (+ 1 PhD student and 1 master student)
- R&D: Development of new resins, techniques and applications
- Products used in several domains

Radiopharmacy
and
Nuclear Medicine

Environment and
Bioassay

Geochemistry
and
Metals Separation

Decommissioning




Research interests - Radiopharmacy

• Radionuclide production/purification

- Resin and method development 'cold'
 - Cooperation with cyclotrons & reactors (NL, RN producers,...)
 - Equipment provider (targetry, synthesizer,...)
- Separation of radionuclides from irradiated targets
 - Diagnostics: Zr-89, Cu-61/4, Ga-68, Ge-68, Ti-44/5, Tc-99m, Sc-43/4...
 - Therapy: Ac-225, Lu-177, Tb-161, Cu-67, Pb-212, Sn-117m, Sc-47...
- Requirements for resins:
 - No selectivity for target material, high selectivity for product
 - Elution under 'soft' conditions in small volume => labelling/injection
 - Fast kinetics
 - Combining several resins can facilitate the separation
 - Conversion (high acid to dilute acid)
 - Removal of impurities upfront

Radiopharmacy
and
Nuclear Medicine





Research interests - Radiopharmacy

Radiopharmacy
and
Nuclear Medicine

- Quality control
 - Cartridge based methods (e.g. Sr-90 in Y-90,...)
 - New option “TK-**El**Scint cartridges” (impregnated plastic scintillator beads)
 - “**Sheets**”
 - p.ex. DGA sheets (functionalized TLC for Ra-223, Ga-68, Pb-212,... => CVUT Prague), CU iSheets,...
- Decontamination of effluents/waste (Ge-68, lanthanides, radioiodine,...)
- Radiolysis stability (polymer, radical scavengers,...)
- Determination of radionuclides (mainly used in therapy, generally Lu-177 and Ac-225) in environmental and bioassay samples



Cu-61/4 separation on TK201

- Cu-61/4 separation from solid Ni targets
 - Target dissolution in high HCl => 6M HCl
 - TK201 retains Cu, Zn, Co, Fe, Ga at 6M HCl
 - Difficult to get a clean Cu fraction in dilute HCl
 - Run through TBP (or TK400) for Fe/Ga removal
 - Separation of remaining radionuclides on TK201
 - Preferably avoid water (risk of Zn co-elution)

Svedjehed et al. *EJNMMI Radiopharmacy and Chemistry* (2020) 5:21
<https://doi.org/10.1186/s41181-020-00108-7>


EJNMMI Radiopharmacy
and Chemistry

RESEARCH ARTICLE

Open Access


Automated, cassette-based isolation and formulation of high-purity [⁶¹Cu]CuCl₂ from solid Ni targets

Johan Svedjehed¹, Christopher J. Kutryk², Jonathan W. Engle^{2,3} and Katherine Gagnon^{1*}



- Svedjehed et al. use of NaCl/HCl for better pH control of eluate
- Also being used for Zn separation
- Not applicable to solid Zn targets


Cu Purification on TK201 after sur TBP for Fe & Ga removal






Upcoming: TK250 Resin

- CU Resin high selectivity for Cu over Zn but loading from pH >2 required
- Difficult to automatize in case of solid Zn targets
- Upcoming TK250 Resin:
- Cu retention from low acid up to 3M HCl
- No selectivity for Ni **and** Zn
 - Tested up to 300mg each
- Cu elution in 6M HCl
- Rather low Cu capacity (~0.13mg/g)




Separation on 2 mL TK250 Resin



Cu separation from 300 mg Ni on 2mL TK250 Resin


Separation on 2 mL TK250 Resin




Cu separation from 300 mg Zn on 2mL TK250 Resin



Upcoming: TK250 Resin



- D_f typically $>10^3 - 10^4$
- 6M HCl to low HCl on TK201 Resin
- Next steps:
 - Optimisation of resin composition
 - Upscale and stability testing
- Integration in sequential separation scheme for Cu and Ga from Zn targets e.g. with TBP Resin (for Ga)





Zr-89 chloride via TBP and TK400

Nuclear Medicine and Biology 136–137 (2024) 108943



ELSEVIER

Contents lists available at ScienceDirect

Nuclear Medicine and Biology

journal homepage: www.elsevier.com/locate/nucmedbio



$[^{89}\text{Zr}]\text{ZrCl}_4$ for direct radiolabeling of DOTA-based precursors[☆]

Serge K. Lyashchenko^{a,b,*}, Tuan Tran^a, Steffen Happel^c, Hijin Park^a, David Bauer^b, Kali Jones^b, Tullio V. Esposito^b, NagaVaraKishore Pillarsetty^b, Jason S. Lewis^{a,b,d}

^a Radiochemistry and Molecular Imaging Probe Core Facility, Memorial Sloan Kettering Cancer Center, New York, NY, USA

^b Department of Radiology, Memorial Sloan Kettering Cancer Center, New York, NY, USA

^c Trikem International, Inc., USA

^d Program in Molecular Pharmacology, Memorial Sloan Kettering Cancer Center, New York, NY, USA




Table 1. Summary of Measured Iron Content in TBP-purified solutions.

Purification Intervention	Measured Iron Content (ppm)	Source
No TK400, TBP only	32.7–38.8 (n = 6)	Graves et al.
Single TK400, followed by TBP	8 (n = 3)	This Study
Double TK400, followed TBP	< 1 (n = 3)	This study

- Improvement of method published by Graves et al. (TBP only) => insufficient Fe removal
- Load and rinse on TBP Resin at ~10M HCl, elution in dilute HCl
- Use of 2xTK400 before TBP Resin for Fe removal
- Production of 11.1 – 14.4 GBq of $[^{89}\text{Zr}]\text{Zr-PSMA-617}$ and $[^{89}\text{Zr}]\text{Zr-PSMA-I&T}$
- Apparent specific activities of 11.1 - 14.4 MBq/ μg
 - 2–3x more than before at industrial quantities.
- On-going:
 - Use of TK201 instead of TK400 for impurities removal (catch additional impurities? e.g. Cu)
 - Alternative methods for Zr oxalate conversion to Zr chloride (avoiding QMA)

Table 2. Summary of Radionuclide Purity Measurements in $[^{89}\text{Zr}]\text{ZrCl}_4$ Solution.

Batch	$[^{89}\text{Zr}]\text{ZrCl}_4$ - Batch 1	$[^{89}\text{Zr}]\text{ZrCl}_4$ - Batch 2	$[^{89}\text{Zr}]\text{ZrCl}_4$ - Batch 3	$[^{89}\text{Zr}]\text{ZrCl}_4$ - Batch 4
Radionuclidic Purity	$\geq 99.9\%$	$\geq 99.9\%$	$\geq 99.9\%$	$\geq 99.9\%$
% of ^{88}Zr	$6.9 \times 10^{-10}\%$	$2.9 \times 10^{-10}\%$	$4.7 \times 10^{-9}\%$	$1.2 \times 10^{-8}\%$
% of ^{88}Y	$3.6 \times 10^{-10}\%$	$2.0 \times 10^{-10}\%$	$2.2 \times 10^{-9}\%$	$5.1 \times 10^{-9}\%$




Ga separation on TK400 Resin

TK400 Resin => use for Ga solid Zn targets (alternative to ZR Resin)

- Ga retention on TK400 from high HCl, elution in low HCl
- No Zn retention
- Faster kinetics than ZR Resin

W. Tieu et al. use of single TK400 cartridge for solid Zn targets (purity lower than ZR/TK200)

Svedjehed et al. use of TK400/A8/TK200 (all 1 mL) for solid Zn targets



On-going:

- Use of TK201 instead of TK200 (idea by Bryce Nelson)
 - Potentially better Zn removal
- Improvement of TK400 Resin / Ga yields (Svedjehed: ~80%)

Demystifying solid targets: Simple and rapid distribution-scale production of [^{68}Ga]GaCl₃ and [^{68}Ga]Ga-PSMA-11

Johan Svedjehed, Martin Pärnaste, Katherine Gagnon*

Cyclotrons and TRACERcenter, GEMS PET Systems AB, GE Healthcare, Uppsala, Sweden



At-211 separation

- Requests for cartridge based separation of At from Bi targets in HNO_3 . Resin approach already used by Burns et al. (3-octanone)
- Eriksen et al. showed At separation from Bi possible in HNO_3 via LLX using Octanol (\Rightarrow TK400 Resin)
- Tereshatov et al. tested several extraction chromatographic resins for At separation from Bi incl. TK400
- At elution via alcohol (removal of org. phase + At from resin)
- TK400 and three additional resins currently being tested (“TK401”, “TK402”, TK200) – standard and new support
- Elution via NaOH possible?
- Currently shipping samples
- Also working on resins for Rn-211/At-211 generator



Chemical Engineering Journal
Volume 464, 15 May 2023, 142742



Mechanism of astatine and bismuth sorption on extraction chromatography resins from nitric acid media



Under development: range of impregnated membrane filters

Upcoming: impregnated membrane filters ('Discs')

Several Discs under development and beta testing:

➤ TK100, TK201, CU , GA, TK101...

➤ 25mm and 47mm diameter

First Disc in range: **TK-GrossAlpha Disc**

=> Alpha measurements

pH 1 - 2, typically up to 100mL samples



Filtration (1 – 10mL/min)

High retention of alpa emitters

- Mainly retained on the surface
- Suprisingly good resolution

Discs glued on steel support => alpha spec

Presence of Thor Ra in QC samples?





DGA Sheets/iSheets




TO-DGA (normal DGA) and TEH-DGA (branched DGA) impregnated TLC paper

- Developed at CVUT (Kozempel et al.)


QC of radionuclides and generator eluents (p.ex. Ra-223, Ac-225/Bi-213, Pb-212, Ge-68/Ga-68 ...)

- TLC scanner or radiometer/LSC or HPGe after cutting


Run under acidic conditions => radionuclidic purity



A scheme of chromatographic separation of mixture of ^{227}Ac and his daughter's nuclides. ^{227}Th remains on start, ^{227}Ac has the retention factor ca 0.2, ^{211}Pb ca 0.7 and ^{223}Ra ca 0.9.



Radiochromatogram measured immediately after separation. Low abundant radionuclides of ^{227}Ac were not detected.



Radiochromatogram measured one hour after separation. Decay and ingrowth of ^{211}Pb is clearly visible.

- Now also available based on iTLC support (faster development, higher DGA load)

=> **DGA iSheets (iSheets = based on iTLC paper)**

- 2D TLC for radionuclide screening?



CU iSheets

Poster presented at Terachem 2022

(Svedjehed et al.)

QC of Cu radiolabeled peptides (labeled vs free Cu)


- Shown: [^{61}Cu]Cu-NOTA-octreotide

Spotting/run on three different papers after labeling:

- Whatman and iTLC without modification and
- CU extractant impregnated iTLC paper.

Both iTLC paper (impregnated/non-impregnated) developed in less than 10min, Whatman took 25 – 30 min.

CU extractant impregnated iTLC paper showed superior resolution



- Other systems under development /testing
- Next: TK213 (Ac, Lu)



Some other on-going projects

- Ac purification via TK221 (or TK222)
- Ra purification and recycling
 - TK101/TK102 (Ra/Ba separation)
 - Development of new Ra Resins
- Upscale of radiolanthanide separations (multi-gram):
Lu-177 and Tb-161
- Tc-99 via cyclotron (TK202, C8 and AlOxA)
- Other radiometals
 - **Ge, Hg, Ag, Pd,...**
- Rapid QC
 - Discs, TK-**El**Scint cartridges (plastic scintillator beads based)
- Decontamination
 - Effluents and reaction wastes
 - I, Ge, Lu, O-18...
- Fate' of RN in the environment
 - Separation methods
 - Mainly longer lived RN (=> therapy)
 - Ac-225/7, Lu-177(m), radioiodine,...
 - Quantification
- Rn-211/At-211 generator
- Microfluidics
- Other 'geometries' & 'Non-resin' separation materials
- Hydrometallurgy
- Analytical applications

Thank you for your attention!

Interested in collaborations?
shappel@triskem.fr

SUBSCRIBE TO OUR NEWSLETTER

To keep updated with our latest developments, news and agenda for a year, subscribe to the TrisKem Info here





TK202 Resin


Tc retention from high NaOH (5 - 7M)

- Dissolved Mo targets
 - Increased Tc (Re) retention at higher Mo concentration
- Clean separation from other elements tested


Re used as homologue

Elution in small volume of water

- Eluate will still alkaline and will contain Na
- Pass through CEX for 'neutralisation' and Na^+ removal and through aluminium oxide for trace Mo removal and recovery as 0.9% NaCl solution



Re/Tc separation from Mo on TK202 Resin




Re separation from selected elements on 2 mL TK202 Resin cartridge, load and rinse at 1 BV/min, elution at 0.25 BV/min



Tc-99m via cyclotron route

TK202 (2mL column) - Mo/Re separation - 2g/2 μ g - load from 5M NaOH




Tests performed cold with 2g Mo and 2 μ g Re


- 2 mL TK202 cartridge
- 2 mL C8 cartridge
- 1 mL AlOx cartridge

Method similar to Zeisler et al.
High Re yield (~90%) in 2 – 3 mL 0.9% NaCl solution

Tc fraction acidification and Na removal on 2mL C8 cartridge




Trace Mo removal on Al₂O₃ cartridge (1ml cartridge)






Tc-99m from large Mo targets




Test with 200g Mo
 ~160 mL TK202 column
 Load from 6 - 7M NaOH - elution in water
 Pass through C8 cartridge for acidification and Na removal
 Final concentration/conversion to 0.9% NaCl or 8 mL AlOxA cartridge





Tc-99m separation from Mo targets – suggested scheme (similar to Zeisler et al.)



- ① Pre-cond. TK202 – 5-7M NaOH → alkaline waste
- ② Load Mo/Tc on TK202 → Mo recovery
- ③ Rinse TK202 – 5-7M NaOH → Mo recovery
- ④ Purge TK202 – Air → Mo recovery
- ⑤ Pre-cond. SCX – HCl then H₂O → Aq. waste
- ⑥ Pre-cond. Al₂O₃ – H₂O → Aq. waste
- ⑦ Elute Tc from TK202 on SCX and load on Al₂O₃ – H₂O
- ⑧ Elute Tc from Al₂O₃ – NaCl 0,9% → Tc recovery

TK202 : 35-75 or 75-150μm
 X : 3-ways valve
 X : 4-ways valve
 SCX : Strong Cation Exchange
 Al₂O₃ : Acidic Alumina

Developed with ReO₄⁻ as TcO₄⁻ surrogate


Re recovered on saline solution from alkaline

Separation with 2g Mo → From 20mL to 2mL

Separation with 200g Mo → From 3L to 20mL



Ag-111 from Pd targets – LN/TK200



Tosato et al.
EJNMMI Radiopharmacy and Chemistry (2023) 8:43
<https://doi.org/10.1186/s41181-023-00232-0>

EJNMMI Radiopharmacy
and Chemistry

RESEARCH ARTICLE

Open Access



Chromatographic separation of silver-111 from neutron-irradiated palladium target: toward direct labeling of radiotracers

Marianna Tosato^{1,2}, Andrea Gandini³, Steffen Happel⁴, Marine Bas⁴, Antonietta Donzella^{5,6}, Aldo Zenoni^{5,6}, Andrea Salvini³, Alberto Andrichetto⁷, Valerio Di Marco² and Mattia Asti^{1*}