SOMMAIRE

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• Bibliographie

- (1) Coma et al. "Automated separation of ⁹⁹Tc using plastic scintillation resin PSresin and openview automated modular separation system (OPENVIEW-AMSS)", Journal of Radioanalytical and Nuclear Chemistry (2019) 321:1057–1065. <https://doi.org/10.1007/s10967-019-06659-7>
- (2) Barrera et al. "A new plastic scintillation resin for single-step separation, concentration and measurement of technetium-99" Analytica Chimica Acta 936 (2016) 259-266. <https://doi.org/10.1016/j.aca.2016.07.008>
- (3) Tarancón et al. "A new plastic scintillation resin for single-step separation, concentration and measurement of ⁹⁹Tc", presented at the NRC9 (29/08/16 – 2/09/16, Helsinki, Finland)
- (4) Hidex eBook "Liquid Scintillation Measuring Procedures: New Developments" <https://hidex.com/ebooks/liquid-scintillationmeasuring-procedures/measuringprocedures/radionuclides-from-nuclearfission-activities/2-3-14-tc-by-rad-disk-andpsresin/>
- (5) J. Garcia & A. Tarancon, "Radionuclide determinations with PS Resin MASS WaterRadd", presented at the European Users Group Meeting in Cambridge (UK) - 21/09/2018, https://www.triskeminternational.com/scripts/files/5bae2550c30ed4.50583030/11_j-garcia_atarancon_radionuclide-determinations-withps-resin_mass_waterradd.pdf
- (6) Giménez et al. "A new method based on selective fluorescent polymers (PSresin) for 3 the analysis of ⁹⁰Sr in presence of ²¹⁰Pb in environmental samples" Applied Radiation and Isotopes, 199, 2023, <https://doi.org/10.1016/j.apradiso.2023.110879>
- (7) Tarancón et al. "Recent applications of Plastic Scintillation Resins », presented at Raddec-Triskem Workshop 2024, (18/04/24, Portsmouth, UK)

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- 24th International Conference on Radionuclide Metrology and its applications (ICRM 2025), Paris (France), 19/05 – 23/05/2025 <https://icrm2025.org/>

- Congreso SEMNIM XLI, Valladolid (Espagne), 21/05 – 23/05/2025 <https://semnim.es/congreso/congreso-semnim-xli-valladolid/>

- 55th Annual Scientific Meeting of the Australian and New Zealand Society of Nuclear Medicine (ANZSNM), 23/05 – 25/05/2025, Melbourne (Australie) <https://www.anzsnm.org.au/eventdetails/27648/anzsnm-55th-annual-scientific-meeting>

- 11. RCA-Workshop, 03/06 – 05/06/2025, Dresden-Rossendorf (Allemagne) <https://www.vkta.de/veranstaltungen/11-rca-workshop-vom-03-05-juni-2025-in-dresden-rossendorf/>

- Canadian Radiotheranostics Leaders' Summit 2025, 12/06 – 13/06.2025, Toronto (Canada) <https://www.canadianisotopes.ca/2025-leaders-summit/>

- 2025 SNMMI Annual Meeting, 21/06 – 24/06/2025, New Orleans (USA) <https://snmmi.org/AM/AM/Home.aspx>

- Conference Goldschmidt 2025, 06/07 – 11/07/2025, Prag (Czech Republic) <https://conf.goldschmidt.info/goldschmidt/2025/meetingapp.cgi>

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
Figure 1. TLC scans of $^{61}\text{Cu}\text{Cu-NOTA-octreotide}$ spotted on: top, iTLC-SG; middle, Whatman paper; bottom extractant-impregnated iTLC-SG. Example 1 notes elevated levels of unlabelled ^{61}Cu , while example 2 notes comparable levels of labelled to unlabelled ^{61}Cu .

Figure 1: Scans TLC de $^{61}\text{Cu}\text{Cu-NOTA-octreotide}$ spottés sur: en haut, iTLC-SG; au milieu, papier Whatman; en bas, CU iSheets.

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• Bibliographie

[1] J. Svedjehed et al.: "New extractant impregnated iTLC-SG paper facilitates improved TLC analysis for Cu radiolabelled peptides", poster presented at TERACHEM 2022, 14 – 17 September 2022, Bressanone (Italy). Poster available online on our website.

[2] J. Svedjehed, K. Gagnon. A quest for simplicity: Automated cassette-based purification of $^{61}\text{Cu}\text{CuCl}_2$ from solid Ni targets using a single time-list. Nucl Med Biol, 108-109, S1 (2022), P-220, ppS170.

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
Coming soon : Extractive Discs

- Range of selective impregnated membrane filters
- Ø = 25mm and 47mm Discs will be available
- Various extractant systems (TK100, TK201, GrossAlpha, CU,...) in beta-testing
- More Discs under development
- For use in filtration or passive sampling set-ups

Use of TK-GrossAlpha Discs for direct alpha spectrometry of aqueous samples

Alpha spectrum, 25mm TK-GrossAlpha Disc,
20 mL HNO₃ at pH2 spiked with
30 mBq each Am-241 and Pu-239

Disc preparation

Sample filtration
(at pH 1 – 2)

↓


Drying of the disc

↓

Glueing on a
rigid support


↓

Alpha spectrometry
measurement

Alpha Spectrum,
47mm TK-GrossAlpha Disc,
100 mL Pedras mineral water at pH 1

TK100 Discs for passive sampling in soils


DGT set-up using TK100 Discs for Sr and Pb capture

« Selective Diffusive Gradients in Thin Films (DGT)
for the Simultaneous Assessment of Labile Sr and Pb
Concentrations and Isotope Ratios in Soils »

Anal. Chem. 2022, 94, 6338–6346

Stefan Wagner, Jakob Santner, Johanna Irrgeher, Markus Puschenreiter,
Steffen Happel, and Thomas Prohaska

TK201 Discs for Tc-99 preconcentration from waste water

Experimental setup, Tc-99 concentration and purification from
waste water using TK201 discs for preconcentration

« Quantification of technetium-99 in wastewater by
means of automated on-line extraction chromatogra-
phy – anion-exchange chromatography – inductively
coupled plasma-mass spectrometry »

J. Anal. At. Spectrom., 2024, 39, 2774-2782

Maximilian Horstmann, C. Derrick Quarles Jr, Steffen Happel,
Michael Sperling, Andreas Faust, David Clases* and Uwe Karst