

The background of the slide features a close-up, slightly blurred image of several laboratory glassware items, including graduated cylinders and vials, arranged diagonally. The lighting is warm, with a yellowish-orange glow on the left side, creating a professional and scientific atmosphere.

Overview and new Developments Radioanalytics

RRMC Workshop 2025

Steffen Happel

03/11/2025



TrisKem International



- Based in Rennes (France)
- Independent company since 2007 (before part of Eichrom Europe)
- Main product line: extraction chromatographic resins and other separation materials
- Staff : 25
- R&D and TechSupport group:
 - 3 RadChem PhD, 2 Technicians (+ 1 PhD student)
- R&D: Development of new resins, techniques and applications

Radiopharmacy
and
Nuclear Medicine

Environment and
Bioassay

Geochemistry
and
Metals Separation

Decommissioning




Cartridge packing

1mL and 2mL cartridges packed fully automatically

Print/label on cartridge body replaced by laser print on cap

Before:



After:





TK Resin naming

There is some reason...

Main domain of application

- TK1NN: Environmental monitoring, bioassay,..
- TK2NN: Radiopharmacy
- TK3NN: Analysis of decommissioning samples and nuclear waste
- TK4NN: Geochemistry
- TK5NN: Hydrometallurgy
- TK6NN: 'custom resins'

Functional groups:

- TK1**00**/1/2 => Crown-ether
- TK2**21**/2/5/7 => DGA
- TK2**11**/2/3 => « organophosphor » acids



TK200 Resin

Based on TOPO extractant


High retention of actinides

Applications:

- Use for very efficient U removal from Pu
- Determination of Tc-99 in water samples
- Ga-68 production (in combination with ZR Resin or TK400 Resin)
- Actinide separation from water samples

Extracts actinides even at pH 1 - 2 (nitric acid)


- Preconcentration and purification of selected actinides on same column
- 'In the field'?






Actinides on TK200

(all data Van Es et al.)




vanEs et al. $[HNO_3]$ M



VanEs et al. $[HCl]$ M



U/Th separation on TK200



Load: 3M HNO₃ or ≥ 1 L pH2 (HNO₃)

Very clean U/Th separation

Oxalate instead of carbonate



TK200 Resin - U/Pu separation

Recent publication by Huang et al.

Better U removal: $D_f(U) > 10^9$

Additional U removal via $He+NH_3$

Overall $D_f(U) > 10^{13}$

Pu isotopes incl. Pu-238 via ICP-MS/MS





Figure S1. Analytical procedure for determination of plutonium isotopes (^{238}Pu , ^{239}Pu , ^{240}Pu , ^{241}Pu)




Actinides on TK200 – Application


Elution of Th from 1g TK200 resin cartridge




Elution of Pu from 1g TK200 resin cartridge



Elution of U from 1g TK200 resin cartridge



- Data by N. Vajda (RadAnal)
- Method development on-going
- Direct load of U, Th and Pu from acidified water samples (here 900 mL on 2 mL TK200)
- Sequential separation on TK200
- Automatisations & 'in the field' preconcentration'?





TK200 – direct Pu load/separation

One TK200 cartridge

- Preconcentration and purification

Automized separation

Acidified water samples (1 L)

Flow rate 15 mL/min

$D_F(U): 10^4 - 10^5$

LoD:

- 0.32 $\mu\text{Bq/L}$ Pu-239
- 2.00 $\mu\text{Bq/L}$ Pu-240




Talanta

Volume 262, 1 September 2023, 124710



A novel strategy for Pu determination in water samples by automated separation in combination with direct ICP-MS/MS measurement

Yuyi Ni^a, Wenting Bu^a, Ke Xiong^a, Sheng Hu^a, Chuting Yang^a, Liguao Cao^b







TK221 Resin

Resin based on a mixture of diglycolamide and phosphine oxide + traces long chained alcohol on inert support with aromatic groups.

- Main applications in RadPharm:
 - Lu, Tb concentration from high acid and elution in small volume of dilute HCl
 - Ac-225 purification
 - Improved radiolysis stability
- Potential use in actinide separations?



(Papp, I., Vajda, N. & Happel, S. An improved rapid method for the determination of actinides in water. *J Radioanal Nucl Chem* **331**, 3835–3846 (2022). <https://doi.org/10.1007/s10967-022-08389-9>)




TK221 Resin – actinide separation

- Synergies between extractants
- Higher U retention than DGA
- Higher Am retention than TRU

Cooperation with Nora Vajda

Method development for water samples


- TEVA/TK221 method
- Ideally later also soil and decommissioning samples
- Influence of Fe and Ca?
 - Ca leads to stronger Am retention



Papp, I., Vajda, N. & Happel, S. *J Radioanal Nucl Chem* (2022). <https://doi.org/10.1007/s10967-022-08389-9>



TK221 Resin



Th and Pu removed with TEVA.

Am and U separation on TK221 Am elution before U

Table 3 Recovery of actinide tracers from spiked water samples

	Actinides determination	
	Without Np separation	With Np separation
	Yield	Yield
	%	%
TAP water		
²³⁰ Th	90 ± 8	86 ± 7
²³⁹ Pu	108 ± 7	95 ± 7
²³⁷ Np	—	91 ± 9
²⁴¹ Am	103 ± 7	97 ± 6
²³³ U	103 ± 7	70 ± 7
SEA water		
²³⁰ Th	71 ± 7	61 ± 6
²³⁹ Pu	91 ± 7	87 ± 6
²³⁷ Np	—	93 ± 8
²⁴¹ Am	89 ± 7	92 ± 6
²³³ U	88 ± 7	59 ± 6

- Method tested on spiked tap and sea water samples
- High yields (88+ for U and Am)
- Analysis of IAEA-TEL-2021-03 WWOPT successful
- On-going: use for solid samples
- TK201/TK221, TK200/TK221 methods?


Analyte	Target values		Measured values			Relative bias	MARB ^a	Z-score ^b	Test evaluation
	Mean activity concentration	Standard deviation (sd)	Activity concentration	Standard uncertainty	Relative standard uncertainty				
	Bq/kg	Bq/kg	Bq/kg	Bq/kg	%				
²³⁹ Pu	5.93	2.27	5.09	0.24	4.7	14	25	0.37	Accepted
²⁴¹ Am	4.85	0.57	4.73	0.15	3.2	2.5	30	0.21	Accepted
²⁴⁴ Cm	7.02	2	7.19	0.34	4.7	2.4	25	0.09	Accepted

^aMaximum Acceptable Relative Bias

^bZ = |X_{reported} - X_{target}| / s_{dtarget}



Tandem TK200/TK221



- Ling Zhang, Emilia Vassileva, Determination of ultra-trace level ^{241}Am in marine sediment and seawater by combining TK200-TK221 tandem-column extraction chromatography and SF ICP-MS, Talanta, 271, 2024, 125724, <https://doi.org/10.1016/j.talanta.2024.125724>

Conclusion:

For DGA separation, the decontamination factors (DFs) for Pu, U, and Th were calculated to be $(1.8 \pm 0.4) \times 10^3$, $(3.0 \pm 1.0) \times 10^4$, and 24 ± 6 , respectively.

For TK221 separation, the DFs for Pu, U, and Th were $(1.5 \pm 0.2) \times 10^4$, $(2.1 \pm 0.3) \times 10^4$, and $(1.2 \pm 0.1) \times 10^3$, respectively.


It can be seen that the TK221 resin displayed remarkably better performance for the removal of Pu and Th than DGA resin, while with an excellent decontamination ability of U, close to DGA.

=> Use for Ra-226/8 instead of DGA?





Ac separation => RP presentation

Analytical questions: bioassay (e.g. BfS) and environment



Load from 2 – 6M HNO₃
 Rinse with 6M HNO₃
 10M HCl => Bi removal and
 0.05M HNO₃ (Fe removal)
 Ac elution in 0.05M HCl
 Po remains on TK221





Sr-90 in seawater – KAERI method

LSC 2024
Raddec/Trsikem Jointworkshop, 18th Apr 2024, Porstmouth

A simple and straightforward technique for analyzing radionuclides in seawater



18 Apr, 2024

Hyuncheol Kim (hckim3@kaeri.re.kr), Gahyun Kim (ghkim97@kaeri.re.kr)



한국원자력연구원
Korea Atomic Energy Research Institute



Sr-90 in seawater – KAERI method

Materials and Apparatus



AMP-PAN (or KNiFC-PAN)

DGA resin

2 mL column



SALT-100 (*WITHTECH Ltd. ; South Korea*)

Eight peristaltic pump

Flow rate: 10 – 100 mL min⁻¹

Applicable with 2 mL / 5 mL column

https://www.withtech.co.kr/en/busi/new_busiList_5.php



Hidex **Q-ARE**

Automated Radionuclide Extraction System

The most advanced automated radionuclide extraction chromatography system dedicated to radionuclide separation from environmental, food and decommissioning samples.

Quick and easy-to-use unattended radionuclide extraction

User friendly, intuitive and hassle free.



<https://www.hidex.com/>

Q uick
A utomated
R adionuclide
E xtraction



Sr-90 in seawater – KAERI method

⁹⁰Sr in seawater; procedure



- Separation & sample prep ~4h (w/o measurement)
- Typically >85% yield
- Presentation by Hyuncheol Kim at the TrisKem / Raddec Workshop, 18.04.2024
- Presentation and video available on our website



Sr-90 in seawater – KAERI method

- Rapid automatized separation of Y-90 from sea water (60 kg)
 - Aim: determination of Sr-90
 - Faster than standard methods
- Potential drawback: Sample needs to be adjusted to 3M HNO_3
=> large amounts of conc. HNO_3
- Ongoing: TO-DGA based resin to allow loading from lower HNO_3
=> 'TK227 Resin'
- Similar approach for Am possible?
- SALT / Q-ARE approach also for used Cs-134/7 via AMP-PAN





TK227 Resin

Addition of ionic liquid increased Y retention (idem e.g. lanthanides, especially heavy lanthanides) at lower acid concentrations, especially in HNO_3 . Also Includes phase stabilizer.


Remark: too high amounts of IL will lead to low elution yields ($\Rightarrow$ TK225 – large amount of IL – use in decontamination of effluents from radiolanthanides).

Remark: Separation on TK227 optimized for mg amounts of carrier.







TK227 Resin – drinking water



Elution study, 1L tap water, loading from 0.1M HNO_3 , Y separation, 2mL TK227 Resin cartridge



Elution study, 1L tap water, loading from 0.1M HNO_3 , Y separation, 2mL DGA,N Resin cartridge



Elution study, 1L tap water, Loading from 0.1M HNO_3 , Y separation, 2mL TK227 Resin cartridge, selected elements

Y retention on TK227 from 1L 0.1M HNO_3

Loading at 10mL/min


Separation from interferences following Kim et al.

Rapid method for Sr-90 in water?


Use of scintillating beads (PSm)?



TK227 Resin – sea water



Elution study, 1L sea water, loading from 0.2M HNO₃, Y separation, 2mL TK227 Resin cartridge



Elution study, 1L sea water, loading from 0.2M HNO₃, Y separation, 2mL DGA,N Resin cartridge




- Sea water requires high HNO₃ concentration for Y retention (at least 0.2M, better higher)
- Loading at 10mL/min
- On TK227 Resin Y breakthrough when loading 1L (<20%) but less than DGA,N (~70%) => higher HNO₃ required to improve yields
- TK227 Resin use in KAERI method showing better yields but acidification to 2M HNO₃ still needed.



TK100/1 Resins

Based on same crownether as SR Resin

- Different solvents, Sr and Pb uptake also between pH 2 and 7
- Concentration and purification on same column




Typical applications:

- Pb-210 in water samples (up to > 5L per 2 mL column/cartridge)
- Sr-90 by ICP-MS (very high Zr-90 decontamination) => NPL
- Ra-226 by ICP-MS => NPL (Load and purification in one step)
 - Agilent application note



TK100 Resin

Development for measurement of ^{90}Sr and ^{226}Ra by Russel and Van Es from NPL



Rapid Analysis of Radium-226 in Water Samples by ICP-QQQ

Application Note
Nuclear, environmental

https://www.agilent.com/cs/library/applications/8800_ICP-MS_5991-8324EN_radium_analysis.pdf

Applied Radiation and Isotopes 126 (2017) 35–39

Authors

Ben Russell¹, Elsje May van Es^{1,2},
Glenn Woods³, David Read^{1,2}

1. National Physical Laboratory,
Teddington, UK

2. Chemistry Department, University
of Surrey, Guildford, Surrey, UK



Contents lists available at ScienceDirect

Applied Radiation and Isotopes

journal homepage: www.elsevier.com/locate/apradiso



Development of an optimised method for analysis of ^{90}Sr in decommissioning wastes by triple quadrupole inductively coupled plasma mass spectrometry

B. Russell*, M. García-Miranda, P. Ivanov

National Physical Laboratory, Hampton Road, Teddington, TW11 0LW, UK



TK100 contains HDEHP



⇒ very high Zr retention ⇒ high Zr decontamination

⇒ Sr elution in $\geq 0.5\text{M HCl}$




TK101 - Radium

TK101 Group 1 and 2



NPL
National Institute of Standards and Technology





- Ra retention from water/dilute acid up to $\sim 0.5M$ HNO_3/HCl
- At higher conc. selectivity closer to SR Resin/TK102 Resin

Data provided by
Russel et al. (NPL)




TK101 Transition Metals



Data provided by
Russel et al. (NPL)
29



TK101 - Ra




Data provided by
Russel et al. (NPL)

- No / extremely low selectivity for Th/U
- Very strong Pb retention => elution in high HCl or citrate



Ra separation on TK101



Good Ra separation when loading from dilute HNO₃/HCl – up to 300mL per 2mL cartridge
When eluting Ra in 3M HNO₃, Ba, Pb, Sr remain retained

Further Ba removal via TK102 possible
Ra/Ba co-eluted in 8M HNO₃
No retention of U, Th, Pt, Ir,...



On-going work: test of other crown ethers and diluents

Aim: two resins: a. improved version of TK101 and b. Ra Resin working at elevated acid concentration and elution in dilute acid

Compound	Name	Abbreviation
	Dicyclohexano-18-crown-6	DCH18C6
	4,4'-(4,5')(5,5')-di-tert-butyl-cyclohexano-18-crown-6	DBDCH18C6
	Dibenzo-21-crown-7	DB21C7
	4,4'-(4,5')(5,5')-di-tert-butyl-dibenzo-21-crown-7	BDB21C7
	Dibenzo-24-crown-8	DB24C8

- Obtained data => Modelisation
- Upcoming: Ra Workgroup with research groups
 - Environment/medical
 - Modelisation
 - Synthesis labs




TK102 Resin

- Modified version of SR Resin
 - Same crown-ether
 - Diluent, inert support and CE/D ratios => different
 - Higher Sr, Pb and Ba retention than SR Resin
 - Less bleeding of organic materials
 - Generally higher D_w and capacity for Sr, Pb and Ba
- Work by Illarion Dohvyi (Poster during ERA14), Marine Bas, Soumaya Khalfallah, Nora Vajda, Steffen Happel
- Optimized for Ra/Ba separation => main use
- Under testing: use in Sr isotope ratios via TIMS (less bleeding) and decommissioning samples




TK102 Resin – Elution curves comparison vs SR Resin - Sr separation




TK102 Resin vs SR resin:
Sr elution study in 8M HNO₃ load
medium

Resins TK102 and SR similar for
the separation of elements
Th/U/Pb/Sr/Ca/Bi/Y/Ca and Ba






TK102 Resin – Elution curves comparison vs SR Resin - Pb separation



TK102 Resin vs SR resin:
Pb elution study with 2M HCl
loading medium


Resins TK102 and SR similar
for the separation of elements
Th/U/Pb/Sr/Ca/Bi/Y/Ca and
Ba






Ra/Ba separation

Separation on 1 mL TK102 Resin (100 - 200 μ m) - ~0.5BV/min



Separation on 1 mL SR Resin (100 - 150 μ m) - 0.5BV/min



- SR Resin: high Ba breakthrough starts after 7 – 8 bed volumes
- TK102 Resin: significantly lower Ba breakthrough
- Suitable for Ba removal from Ra at 3M HNO₃



Tc-99 separation

Tc-99 (difficult to measure – DTM Radionuclide) – 100% beta emitter

TEVA resin allows for Tc separation but quantitative elution needs highly acidic medium

➤ Other options

TK201 Resin

TK202 Resin

TK200 Resin

TK-TcScint




TK201 Resin

Based on tertiary amine (weak Anion Exchanger)
impregnated on inert support



Main applications:

- Tc separation
 - Environmental monitoring
 - Decommissioning
 - Easier to elute
 - Use of NH_4OH or 3M HNO_3
- Cu separation
 - Cu-61/4 from Ni targets
- Zn, Fe separation
 - Geochemistry





TK201 – Actinides





- Re uptake between pH 2 – 5M HCl
- In HNO_3 medium, Re fixed at pH 1-2

- High U and Pu uptake at high HCl
- Elution at low HCl

Data Russel et al. NPL



TK201 – Dw values



Tc well retained at $c(\text{HNO}_3) < 1\text{M}$

Tc retention significantly higher in HCl

No Mo retention at $> 0.7\text{M}$ HNO₃



TK201 Resin – Elution curve



Load of sample at pH 1-2 to retain Re & Tc

Interferences removed during load/rinse

Mo elution at 0.7M HNO₃

Tc/Re elution in 3M HNO₃

Preferable elution option: $\geq 0.5\text{M NH}_4\text{OH}$



Online Solid-Phase Extraction–Inductively Coupled Plasma–Quadrupole Mass Spectrometry with Oxygen Dynamic Reaction for Quantification of Technetium-99


Makoto Matsueda,* Kayo Yanagisawa, Kazuma Koarai, Motoki Terashima, Kenso Fujiwara, Hironobu Abe, Akihiro Kitamura, and Yoshitaka Takagai*



Cite This: *ACS Omega* 2021, 6, 19281–19290



Read Online



JAAS

PAPER


[View Article Online](#)
[View Journal](#) | [View Issue](#)

Check for updates

Quantification of technetium-99 in wastewater by means of automated on-line extraction chromatography – anion-exchange chromatography – inductively coupled plasma-mass spectrometry†

Maximilian Horstmann,^a C. Derrick Quarles, Jr.,^b Steffen Happel,^c Michael Sperling,^d Andreas Faust,^e David Clases^f and Uwe Karst^{*,a}

Cite this: *J. Anal. At. Spectrom.*, 2024, 39, 2774



Second publication including TK201 Discs for higher sample load => better LoD





TK202 Resin

Polyethylene Glycol (PEG) grafted on inert support


Retention of chaotropic anions e.g; TcO_4^- in the presence of kosmotropic anions (SO_4^{2-} , CO_3^{2-} , OH^- , MoO_4^{2-} ,...)

For samples rich in Mo: Tc yield > 90% for 6 – 8g Mo per g TK202

Analytical application: e.g. Tc-99 in concrete samples



D_w values for Tc, Re and Mo on TK202 Resin, at varying NaOH concentrations. Tc data taken from Cieszykowska et al.



D_w values for Tc in 5M NaOH using 40 mg TK202 Resin, increasing amounts of Mo. Data taken from Cieszykowska et al.



TK202 Resin


Retention of Tc from concentrated NaOH medium (5 - 7M)

- Alkaline Fusion e.g. decommissioning samples
- Dissolution of Mo target
- Clean separation from other tested elements


Re can be used as internal standard

Elution in a small volume of water

- Eluat remains slightly alkaline
- Load on CEX to neutralise medium + remove Na^+
- Potentially load on aluminum oxide to remove last Mo traces + elution in 0.9% NaCl



Re/Tc separation from Mo on TK202 Resin




Re separation from selected elements on 2 mL TK202 Resin cartridge, load and rinse at 1 BV/min, elution at 0.25 BV/min.



TK200 Resin

Rapid determination of ^{99}Tc in water samples using $\text{Ti}(\text{OH})_3\text{-TcO}_2$ co-precipitation and TK200 resin by liquid scintillation counting

Ni Yuan^a, Quan An^a, Shan Xing^{a, b}, Xiongxin Dai^a, Xiaolin Hou^{c, d}, Yonggang Yang^a, Yan Ma^a





TK200 Resin

Rapid determination of ^{99}Tc in water samples using $\text{Ti}(\text{OH})_3$ - TcO_2 co-precipitation and TK200 resin by liquid scintillation counting

Ni Yuan^a, Quan An^a, Shan Xing^{a, b}, Xiongxin Dai^a, Xiaolin Hou^{c, d}, Yonggang Yang^a, Yan Ma^a

1L v





Fig. 1. Schematic diagram of the chemical procedure for separating ^{99}Tc from the water sample.

using TDCR-LSC





TK-TcScint



UNIVERSITAT DE
BARCELONA

Plastic scintillating beads impregnated with selective extractant
Developped by University of Barcelona

- García, Tarancón, Bagán
- « TK-ElScint » product line
 - 1st product: « TK-TcScint »
 - Quaternary ammonium + phase modifier (similar selectivity to TEVA)
 - Environment/decommissioning => Tc-99 by LSC
 - New: TK-SrScint
 - Sr and Pb
 - More products under development:
 - gross alpha(/TK-GA), DGA(/TK227), TK101 (Pb, Ra,...),...



Direct measurement of the cartridge by LSC after loading and rinsing

- NO elution/evaporation/aliquoting => easy automatisaton

Chemical yield via Re/ICP-MS in eluates.



TK-TcScint

Use of TK-TcScint in aqueous/urine samples for Tc-99 determination (Garcia et al., TKI UGM Cambridge 2018)

MOP:

2ml cartridge using Vacbox

1mg Re carrier


Precondition with 2ml 0.1M HCl

Load 10ml sample in 0.1M HCl

Rinse 4x2ml DI H₂O

Results

Recovery of Rhenium (by ICP-OES)	> 98.8 %
Recovery of ⁹⁹ Tc (by LS):	> 98.8 %
⁹⁹ Tc Detection Efficiency (%):	89.5(0.6)
Background (cpm):	1.09
Quenching Parameter (SQP(E)):	787(7)





TK-SrScint


Impregnated plastic scintillation beads (PSm)

TK102 extractanten system (better loading volumen, capacity)

Crown ether + fluorinated alcohol

Sr and Pb retention => chemical yield typically Sr or Pb carrier in eluate via ICP

Detection efficiency for Sr-90: >90%, after Y-90 ingrowth >190%



Normalized LSC spectra Sr-90/Y-90 (blue) and Pb-210/Bi-210 (red)
on TK-SrScint at $t = 0$ (A) and after ingrowth of the daughters at $t > 21$ Tagen (B).



TK-SrScint

Well suitable for decommissioning samples (no presence of Pb-210)

Environmental samples, potential for false positives through Pb-210 (and Bi-210)


- No Sr elution, Pb and Sr remain on the resin, Y-90 and Bi-210 ingrowth
- Tarancon et al. Pb-210 removal through co-precipitation with iodate in presence of Ca at ΔT . Sr remains in solution
- Sr co-precipitation with Ca-Phosphate
- Dissolution in 8 M HNO_3
- Load onto TK-SrScint
- Rinse with 8M HNO_3 , followed by 6M HNO_3 and 6M LiNO_3 (acid removal)

For 1L Wasserproben: additional preconcentration via carbonate co precipitation

=> Chemical yield for full process between 63 % and 81 %



TK-SrScint – Rapid method Sr-90 in milk



Sr-90 determination using TK-SrScint



TK400 Resin

Long chained alcohol

- First work by Knight et al. on Np/Pa separation



Retention only at high HCl concentration (>6M HCl), elution in low HCl, water,...

Main application: Pa separation & Ga-68 production


Also retains **Mo**, **Fe**, **Po**, **Nb**,... working on **Sb** and **At**

Higher Fe capacity than e.g. TRU Resin (~15mg Fe/g TK400)

Knight et al.





Data provided by B. Russel (NPL)





Fe/Mo/Nb separation



- Part of sequential method for determination of Fe, Ni, Zr, Mo, Nb, Tc in stainless steel (PhD student with Subatech, cooperation with IPHC)
 - Project includes modified NI Resin for improved selectivity (esp. Cu)
- Recovery of Fe/Nb/Mo from high HCl on TK400
- Majority of other elements removed during load and rinses (9M and 5M HCl)
- Fe/Nb and Mo eluted in dilute HCl => separation on ZR Resin



Calixarene based resins for Cs separation

AMP-PAN and KNiFC-PAN well suitable for Cs concentration from aqueous matrices but:

Cs elution difficult, leading to high matrix Cs containing solutions

- Use of NH_4OH , $\text{Sr}(\text{OH})_2$ followed by AIX and CEX

Use of calixarene based resins instead

Original work: TK300 Resin

- High Cs/Ba selectivity
- Load from water up to 1M HNO_3
- Interference by K
- Low Cs capacity
- Home made calixarene => upscale too difficult



Calixarene based resins

Aim: two resins

- Separation of Cs and Rb from neutral to weak acid and elution with strong acid and vice versa
- Use of ionic liquids or short-chained alcohols,... as diluents

Preparation of >15 test resins (PR) based on commercially available calixarenes

Characterisation of these resins


- D_W values of different element on test resins in HNO_3 and HCl
- Influence of interfering ions (like potassium) on Cs separation
- Breakthrough and full capacities
- Elution tests for Rb and Cs separation

Work performed by Illarion Dovhyi, presented at last TKI/Raddec Workshop (18.04.2024) => visit our website for more data


Decision on which resins will be fabricated end of this year => beta testing⁵⁷




Tested Calixarenes



Calix[4]arene-bis(tert-octylbenzocrown-6, $C_{72}H_{92}O_{12}$ (BOBCalix)




1,3-alt-25,27-Bis(3,7-dimethyloctyl-1-oxy)calix[4]arene-benzocrown-6, $C_{62}H_{82}O_8$ (MAXCalix)




Calix[4]arene-bis[4-(2-ethylhexyl)benzo-crown-6], $C_{72}H_{92}O_{12}$ (BEBHCAlix)



D_W values of selected cations in HNO_3



Acid dependency of D_W for Cs^+ on PR1-13 in HNO_3




Acid dependency of D_W for Rb^+ on PR 1-13 in HNO_3

- Generally, $D_W \text{ Cs} > D_W \text{ Rb}$, other elements (Ag, Al, Ba, Bi, Ce, Co, Cu, Er, Eu, Ga, Hf, La, Lu, Mo, Nb, Nd, Ni, Pb, Re, Sb, Sc, Sn, Sr, Th, U, Y, Yb, Zn, Zr) not retained from HNO_3
- Ionic liquid based resins: very high D_W from 0.01M to $\geq 1\text{M}$ HNO_3
- Other test resins (non IL) low Cs/Rb extraction from low acid, maximum at 2 – 3M HNO_3
- Generally strong decrease of retention at very high HNO_3



Elution tests with selected test resins


Separation on 1 mL
PR8



Separation of Rb and Cs on PR8
(loading in 0,01 M HNO₃)

- Generally high selectivity for Cs and Rb, interferents well removed
- Cs elution required >5 BV 8M HNO₃

Separation on 1 mL PR8



Separation of Rb and Cs on PR8
(loading in 1 M HNO₃ solution)




New/under development impregnated membrane filters

Coming soon : Extractive Discs


- Range of selective impregnated membrane filters
- Ø = 25mm and 47mm Discs will be available
- Various extractant systems (TK100, TK201, GrossAlpha, CU,...) in beta-testing
- More Discs under development
- For use in filtration or passive sampling set-ups



Use of TK-GrossAlpha Discs for direct alpha spectrometry of aqueous samples




TK100 Discs for passive sampling in soils



« Selective Diffusive Gradients in Thin Films (DGT) for the Simultaneous Assessment of Labile Sr and Pb Concentrations and Isotope Ratios in Soils »

Anal. Chem. 2022, 94, 6338-6346
Stefan Wagner, Jakob Sambrer, Johanna Irgeher, Markus Paschenreiter, Steffen Hoppel, and Thomas Prohaska

TK201 Discs for Tc-99 preconcentration from waste water



« Quantification of technetium-99 in wastewater by means of automated on-line extraction chromatography - anion-exchange chromatography - inductively coupled plasma-mass spectrometry »

J. Anal. At. Spectrom., 2024, 39, 2774-2782
Maximilian Horstmann, C. Derrick Charles Jr, Steffen Hoppel, Michael Sperling, Andreas Faust, David Clasen* and Uwe Kuntz



Special membrane filters impregnated with
extractant(s)

25mm and 47mm Discs available

Commercialized:

- TK-GA Discs (Gross Actinide)

Beta Testing:

- TK100 & TK201 (both published), CU,...
- TK200 (U, Th, Pu) und TK221 (U, Th, Am, Pu) => sequential separation (evtl. TK201 (Pu, Th)) => rapid screening?
- Alphaspec: Retention vs. Retention on the surface (70 - >90%)

Under development:

- TK102/TK227 - Rapid methods for Sr, TK101 (Pb, Ra),...




New/under development: impregnated membrane filters

Coming soon : Extractive Discs


- Range of selective impregnated membrane filters
- Ø = 25mm and 47mm Discs will be available
- Various extractant systems (TK100, TK201, GrossAlpha, CU,...) in beta-testing
- More Discs under development
- For use in filtration or passive sampling set-ups



Use of TK-GrossAlpha Discs for direct alpha spectrometry of aqueous samples



TK100 Discs for passive sampling in soils




DGT set-up using TK100 Discs for Sr and Pb capture

« Selective Diffusive Gradients in Thin Films (DGT) for the Simultaneous Assessment of Labile Sr and Pb Concentrations and Isotope Ratios in Soils »

Anal. Chem. 2022, 94, 6338-6346

Stefan Wagner, Jakob Samner, Johanna Irgeher, Markus Paschenreiter, Steffen Hoppell, and Thomas Prohaska

TK201 Discs for Tc-99 preconcentration from waste water



Experimental setup, Tc-99 concentration and purification from waste water using TK201 discs for preconcentration

« Quantification of technetium-99 in wastewater by means of automated on-line extraction chromatography - anion-exchange chromatography - inductively coupled plasma-mass spectrometry »

J. Anal. At. Spectrom., 2024, 39, 2774-2782


Maximilian Horstmann, C. Derrick Quares Jr, Steffen Hoppell, Michael Sperling, Andreas Faust, David Clasen* and Uwe Kuntz

Typical applications of TK-GA Disc:

- pH 1 – 2 (for some actinides up to 3M HNO₃)
- Preconditioning of the disc: 20% EtOH:water, water, then loading medium
- Loading of the sample, flow rate: 1 - 10mL/min,
- Lower flow rates seem to favor better resolution
- 25 mm Discs typically up to 100mL samples,
- 47 mm Discs typically up to 1L (higher volumes possible depending on extractant system)
- Rinse with water and 20% EtOH:water
- Drying and glueing onto steel disc => alpha spectrometry.
- Optional: subsequent Alpha/Beta via LSC



TK-GA Disc – Am/Pu in spiked acidified water




Alpha spectrum, Am-241 & Pu-239, each ~50mBq, 100mL pH2 HNO₃




TK-GA Disc – U/Po in spiked acidified water

→ Solution synthétique simple dopée en ^{238}U et ^{209}Po
(100mL HNO_3 0,01M ; sans ou avec Ca à 0,1g/L ; 30mBq soit 0,3Bq/L)



Résultats : → Eau minérale naturelle : Pedras (100mL)




No Ra Retention in pH1 , Th very well retained. Good resolution



TK-GA Disc – Ra in spiked acidified water

→ Solution synthétique simple HNO_3 (100mL) 0,01M dopée en ^{226}Ra (30mBq soit 0,3Bq/L) (résultats octobre 2024)


Spectre Alpha du GA disc



No Ra retention at Ph 1 and 2. Ra retained at higher pH but no spectrum



TK-GA Disc – Am/Pu in sea water




Alpha spectrum, Am-241 & Pu-239, each ~30mBq, 100mL sea water, pH2 HNO₃

Higher yields for Pu



TK-GA Disc – Am/Pu in 3M HNO₃



Alpha spectrum, Am-241 & Pu-239, each ~50mBq, 100mL 3M HNO₃

High yields for Pu, no Am retention



Next steps

TK-GA Discs:

- Publication upcoming
- Other matrices (e.g. primary coolant) and acid concentrations
- Spectrum deconvolution => Subatech
- More data on other alpha emitters
- Behaviour of Ra
 - Retention at pH3-4 but no alpha spectrum
 - Migration into the disc?

TK221 and TK200 Discs in 3M HNO₃

- Stacking for AC screening? => rapid screening methods
- Will require 6M NaNO₃ rinse before drying for acid removal...
- => looks good cold, alpha spectra?

TK201 and TK101 Discs and methods optimisation

Sr via TK102 and/or TK227 Discs? LSC? GPC?

More discs...



Passive sampling

- Wagner et al.
- Passive sampling via DGT (Diffusive Gradient in Thin films) => 'bio-availability'
- TK100 Discs for Sr, Pb, Zn and Sr-90
- CU Discs for Cu

analytical
chemistry

pubs.acs.org/ac



Article

Selective Diffusive Gradients in Thin Films (DGT) for the Simultaneous Assessment of Labile Sr and Pb Concentrations and Isotope Ratios in Soils


Stefan Wagner, Jakob Santner, Johanna Irrgeher, Markus Puschenreiter, Steffen Happel, and Thomas Prohaska*



Cite This: *Anal. Chem.* 2022, 94, 6338–6346




Read Online





Tc Reconcentration



Experimental setup, Tc-99 concentration and purification from waste water using TK201 discs for preconcentration

- 1L hospital waste water
- Preconcentration on 47mm TK201 Disc
- Elution in 0.5M NH_4OH
- Readjustment of pH to allow loading onto TK201 column (PrepFast)
- Elution in 0.5M NH_4OH concentration/purification step)
- Load onto IC column
- Elution in 0.15M NH_4NO_3 in very small volume
- ~100% chemical yield
- LoD < 1fg/kg

« Quantification of technetium-99 in wastewater by means of automated on-line extraction chromatography – anion-exchange chromatography – inductively coupled plasma-mass spectrometry »

J. Anal. At. Spectrom., 2024, 39, 2774-2782

Maximilian Horstmann, C. Derrick Quarles Jr, Steffen Happel, Michael Sperling, Andreas Faust, David Clases* and Uwe Karst



Production of « industrial » extraction chromatographic resins

Hydrometallurgy, some decontamination work

- Valorisation / critical metals from 'waste'
- Decontamination / valorisation of wastes, effluents or decontamination agents (e.g. acid)
 - Example: Am-241 for isotope batteries
- Number of reuses

Range of different resins being fabricated and tested

Bigger particle size support: ~300 – 800µm

Higher amount of resins requested


- Challenge: supply of raw materials
- Pricing

Increase of production capacity for these resins



Some other on-going projects

- Rapid tests
 - Impregnated PSm resins
 - Range of 'Test sticks'
 - Suitable impregnated support
 - UM => rapide isotope ratio analysis by MS (metallomics): Cu, Ca, Zn,...
 - NPL
 - Uni Barcelona
- Passive sampling (DGT)
- Separation of DTM
 - SE Resin => Ines presentation
 - Zr-93, Fe, Mo, Nb,...
- Fate' of RN in the environment + bioassay
 - Separation methods
 - Mainly longer lived RN (=> therapy)
 - Ac-225/7, Lu-177(m), radioiodine,...
 - Quantification
- In-field preconcentration
 - Impregnated membranes
 - Cartridges
- Microfluidics
- Other 'geometries' & 'Non-resin' separation materials



Thank you for your attention!



SUBSCRIBE TO OUR NEWSLETTER

To keep updated with our latest developments, news and agenda for a year, subscribe to the TrisKem Info here

