

Abstract

Triskem develops extraction chromatographic resins specific for radionuclides and/or group of radionuclides that are used in different fields such as environmental monitoring, decommissioning, nuclear medicine, and more. Depending on the application, the resin format is not always the most efficient option. This is why in the last years Triskem has initiated the development of new formats such as selective discs/filtering membranes and impregnated Thin Layer Chromatography supports. TK100 Discs are used in DGT devices for the accurate quantification of trace-level ($\mu\text{g}\cdot\text{L}^{-1}$) Sr and Pb concentrations and isotope ratios [$\delta_{\text{SRM } 987}^{(87}\text{Sr}/^{86}\text{Sr})$ and $\delta_{\text{SRM } 981}^{(207}\text{Pb}/^{206}\text{Pb})$] in labile, bioavailable element fractions in soils (1,2). TK201 Discs have been used for the preconcentration of Tc-99 from hospital waster waters before purification and quantification by ICP-MS (3). In the same format, TK-GA Discs have been developed to trap actinides from water samples during filtration and thus allowing subsequent alpha spectrometry analysis (4).

TLC supports functionalized via impregnation with selective extractants is another field of development to provide mostly nuclear medicine with fast quality control regarding labelling yields and radionuclidic impurity determination. DGA Sheets have been developed by Kozempel et al. at CTUP for the determination of radionuclidic impurities e.g. in Ra-223, Ac-225 or Pb-212 solution for use in labelling. They further have been used to quantify Ra-Macropa labelling yields. CU Sheets have been developed to facilitate the analysis of Cu radiolabelled peptides.

Further, the development of other selective Discs/filtering membranes and iTLC supports currently being developed are presented here.

Key words: impregnated discs, iTLC, TK100 discs, TK-GA discs, DGA Sheets, CU iSheets.


Extractive/Impregnated discs: Transposition of different resins extractive system onto discs to simplify preconcentration/purification/detection-counting

- Filtration or passive sampling on 25mm or 47mm diameter discs
- Already in-use/published: TK100, TK201 and TK-GA




Typical applications of TK-GA Discs [4]:

- pH 1 – 2 (for some actinides up to 3M HNO_3)
- Disc preconditioning : 20% EtOH:water, water, then loading medium
- Sample Loading, flow rate: 1 - 10mL/min
 - Lower flow rates favors better resolution
 - 25 mm Discs => up to 100mL samples
 - 47 mm Discs => up to 1L (higher volumes possible depending on extractant system)
- Rinsing with water and 20% EtOH:water
- Drying and glueing onto steel disc => alpha spectrometry
- Optional: subsequent Alpha/Beta by LSC

TK100 Discs for passive uptake of Sr and Pb in soils [1][2]





TK201 discs for Tc-99 preconcentration /determination in waste water [3]



Impregnated TLC supports: for Fast Radionuclidic purity QC Checks


DGA Sheets by Kozempel et al. (CTUP) [5][6]

QC of radionuclides and generator eluents (e.g.: Ra-223, Ac-225/Bi-213, Pb-212,...)



CU Sheets by Svedjehed et al. [7]:

QC of Cu radiolabeled peptides (labeled vs free Cu): Both iTLC paper (SG and SG-CU) developed in less than 10 min. while Whatman paper took 25 – 30 min.



On-going studies

- **Extractive membranes** : TK221 and TK200 for Actinides rapid screening // new application using TK200 for Po in Ac-225 // TK101 for Ra and TK102 for Sr // TK227 for Ra, indirect measurement of Sr-90 via Y-90 // TK201 for Tc

Methods developments and optimisation

- **iTLC sheets** : Zr, TK201, TK213,...

2D TLC for radionuclide screening

Bibliography

- [1] S. Wagner, J. Santner, J. Irgeher et al. Selective Diffusive Gradients in Thin Films (DGT) for the Simultaneous Assessment of Labile Sr and Pb Concentrations and Isotope Ratios in Soils. *Anal Chem*, 2022, 94, 16, 6338-6346, <https://doi.org/10.1021/acs.analchem.2c00546>
- [2] S. Wagner, J. Santner, M; Puschenreiter, et al. Prediction of Radiogenic Sr and Pb Isotope Signatures in Plants using Diffusive Gradients in Thin Films. *Analytical and Bioanalytical Chemistry*, 2026, 418, 1541-1555, <https://doi.org/10.1007/s00216-026-06315-6>
- [3] M. Horstmann , C. Derrick Quarles Jr, S. Happel et al. : Quantification of technetium-99 in wastewater by means of automated on-line extraction chromatography – anion-exchange chromatography – inductively coupled plasma-mass spectrometry. *J. Anal. At. Spectrom.*, 2024, 39, 2774-2782, <https://doi.org/10.1039/D4JA00270A>
- [4] C. Bailly et al. Selective impregnated Membrane Filters for Direct Alpha Spectrometric Analysis of Actinides; *Applied Radiation and Isotopes*, in preparation
- [5] A. Bajzicková, L. Kománková, J. Kozempel et al.: Assessment of 223Ra radionuclidic purity by extraction paper chromatography, *CH - Jaderná a kvantová chemie, fotochemie*, 2015, [RIV/68407700:21340/15:00230428](https://doi.org/10.1007/s00230428)
- [6] D. S. Abou, N. A. Thiele, et al.: Towards the stable chelation of radium for biomedical applications with an 18-membered macrocyclic ligand. *Chem. Sci.*, 2021, 12, 3733-3742, <https://doi.org/10.1039/D0SC06867E>
- [7] J. Svedjehed et al.: "New extractant impregnated iTLC-SG paper facilitates improved TLC analysis for Cu radiolabelled peptides", poster presented at TERACHEM 2022, 14 – 17 September 2022, Bressanone (Italy). Poster available on <https://www.triskem-international.com/>