



**TK227 Resin - Uses for preconcentration
and measurement of Sr-90 via Y-90**

Developments of alternative supports

Aude Bombard (abombard@triskem.fr)

Raddec-Triskem Technical Workshop

Portsmouth – 11/09/2026



Summary

- Context of TK227 Resin development
- Characterisation of TK227 Resin
- Applications
- Developments of alternative supports
- Conclusions



Summary

- **Context of TK227 Resin development**
- Characterisation of TK227 Resin
- Applications
- Development of alternative supports
- Conclusions



Context of TK227 Resin development

Measurement of Sr-90 via Y-90 in seawater samples by DGA,N Resin [1][2]

- Y-90 at equilibrium with Sr-90
- No need to wait 14 days for the (re-)measure @equilibrium after Y-90 ingrowth
 - Rapid methods developed to assay ^{90}Sr at very low levels e.g. Method by Maxwell et al.[1] – 1-40L
 - Method by Kim et al.[2] - 1-50L/60kg

[1] Maxwell, S.L., Culligan, B.K., Hutchison, J.B. et al. Rapid determination of ^{90}Sr in seawater samples. *J Radioanal Nucl Chem* **303**, 709–717 (2015).
<https://doi.org/10.1007/s10967-014-3391-8>

[2] Hyuncheol Kim and Gahyun Kim: “A simple and straightforward technique for analyzing radionuclides in sea water” presented at the Raddec/Triskem Workshop (as part of the 2024 LSC conference) on April 18 in Portsmouth (UK). Available online: https://www.triskem-international.com/scripts/files/662fa5023ef284.68970177/raddec_tki-ws-2024_6_kim.pdf



Context of TK227 Resin development

⁹⁰Sr in seawater; procedure





- Kim et al. [2][3][4]
- Gahyun et al. [5]



Elapsed time for analyzing the 60 kg of seawater?

4 hr excepting the measurement time



[2] Hyuncheol Kim and Gahyun Kim: "A simple and straightforward technique for analyzing radionuclides in sea water" presented at the Raddec/Triskem Workshop (as part of the 2024 LSC conference) on April 18 in Portsmouth (UK). Available online: https://www.triskem-international.com/scripts/files/662fa5023ef284.68970177/raddec_tki-ws-2024_6_kim.pdf

[3] Video KAERI Method: "Simple and Straightforward technique for Sr-90 analysis in seawater". Available online: <https://www.youtube.com/watch?v=i2R2b98oTCU>

[4] Hyuncheol Kim, Yoo Gyum Kang, Yong-Jin Lee et al.: "Automated extraction chromatographic radionuclide separation system for analysis of ⁹⁰Sr in sea water", Talanta 217 (2020) 121055. <https://doi.org/10.1016/j.talanta.2020.121055>

[5] Gahyun Kim, Sang-Do Choi, Jong-Myoung Lim, Hyuncheol Kim, Strontium-90 levels in seawater southeast of Jeju Island during 2021–2023, Marine Pollution Bulletin, Volume 193, 2023, 115258. <https://doi.org/10.1016/j.marpolbul.2023.115258>



Context of TK227 Resin development

BUT acidification of 50L seawater @ 3M HNO₃ requires 9,5-10L of cc HNO₃

- Modification of DGA,N Resin to get Y-90 @ lower HNO₃ concentrations

Initial 3M HNO ₃ – 50L	9,5L of cc HNO ₃
1M HNO ₃	3,2L
0,5M HNO ₃	1,6L
0,1M HNO ₃	320mL

- Ideally, 0,1M HNO₃ would be the best result.



Context of TK227 Resin development

How to increase retention of Y (trivalent elements)
@lower HNO_3 concentrations

➤ **use of ionic liquid**

✓ **All work presented about TK227 Resin was performed by I. Dovhyi**



Summary


- *Context of TK227 Resin development*
- **Characterisation of TK227 Resin**
- Applications
- Development of alternative supports
- Conclusions



Characterisation of TK227 Resin

TK227 resin is comprised of


- Inert support
- TO-DGA extractant
- Ionic liquid
- Long chain alcohol






Characterisation of TK227 Resin

K_d values of elements of interest in HNO₃ medium on TK227 Resin



K_d values of elements of interest in HCl medium on TK227 Resin




- Strong uptake of Y on the HNO₃ concentrations tested (same for Eu)
- Uptake of Sr between 1 and 3M HNO₃
- Strong Th uptake above 0,5M HNO₃
- Uptake of U above 2M HNO₃

In HCl,

- No retention below 0,5M HCl
- Retention increase for Y, Eu, Th and U at high HCl concentration



Characterisation of TK227 Resin



- ✓ Developed TK227 allows for higher retention of Y at low HNO₃ concentrations
- ✓ Y-carrier needed (minimum of 1mg).



Summary


- *Context of TK227 Resin development*
- *Characterisation of TK227 Resin*
- **Applications**
- Development of alternative supports
- Conclusions



Applications:


Y breakthrough studies - tapwater

1L samples acidified @ pH 1, Test of Y breakthrough on 2mL cartridges of both TK227 and DGA,N Resins (Fractions measured by ICP-MS)



TK227 Resin:

- Slight Y breakthrough observable starting at 700mL but remain quite low.
- Elution of Y in 10-20mL



DGA,N Resin:


- Y breakthrough starting immediately with loading at 0,1M HNO₃
- Elution of remaining Y in 10-20mL



Applications:


Y breakthrough studies – seawater

1L samples acidified @ 0,2M HNO₃, Test of Y breakthrough on 2mL cartridge of both TK227 and DGA,N Resins (Fractions measured by ICP-MS)



TK227 Resin:

- Y breakthrough observable starting at 300mL. Remain <20%.
- Elution of remaining Y in 10-20mL
- Better Y retention using higher HNO₃ conc. (e.g. 0,5M HNO₃)



DGA,N Resin:


- Y breakthrough starting immediately with loading at 0,2M HNO₃
- Elution of remaining Y in 50mL



Applications:

Elution studies on tapwater sample

1L Tapwater sample acidified @ pH 1, Elution study on 1x2mL cartridge TK227 using Kim et al. Conditions [1] (Fractions measured by ICP-MS)





- Th, U, Sr, Pb, Cs elute during Loading and Rinsing 2
- Quantitative recovery of Y in 40mL 0,1M HCl
- Only one preconcentration/purification step



Applications:

Use for Am separation @ low pH

- Eu (proxi for Am) with high uptake between 0,01 and 0,1M HNO_3



Method development to be tested to preconcentrate/purify Am from low HNO_3 concentration

Load @ 0,1M HNO_3 ; Rinse with 0,1M HNO_3 ; Hypothesis on rinses: *Rinses with 8-9M HCl then 3M HNO_3 /0,3M HF*; Elution Am with 0,1-0,25M HCl

To be used in combination with TK200 Resin for (U,Th, Pu) separation as TK200/TK227 (see Inés Llopart presentation)



Summary

- *Context of TK227 Resin development*
- *Characterisation of TK227 Resin*
- *Applications*
- **Development of alternative supports**
- **Conclusions**



Development of alternative supports: impregnation of different types of membranes

Resins are not always the optimum format of use depending on final application, sample nature, ...

Kinetics of uptake are being quite rapid for most resins

=> Transfer of the impregnation technique from resins to alternative supports

- tests sticks (see presentation of Alex Tribolet) => DTM screening to decide for further chemistry or not
- Filtering membranes,
- Impregnated Thin Layer Chromatography




Development of alternative supports: impregnated discs already in use

Transposition of different resins extractive system onto discs to simplify preconcentration/purification/detection-counting

- Filtration or passive sampling on 25mm or 47mm diameter discs
- Already in-use/published: TK100, TK201 and TK-GA

TK100 Discs for passive uptake of Sr and Pb in soils [1][2]




DGT set-up using TK100 Discs for Sr and Pb capture

« Selective Diffusive Gradients in Thin Films (DGT) for the Simultaneous Assessment of Labile Sr and Pb Concentrations and Isotope Ratios in Soils »

Anal. Chem. 2022, 94, 6338–6346
Stefan Wagner, Jakob Santner, Johanna Irrgeher, Markus Puschenreiter, Steffen Happel, and Thomas Prohaska

TK201 discs for Tc-99 preconcentration /determination in waste water [3]



Experimental setup, Tc-99 concentration and purification from waste water using TK201 discs for preconcentration

« Quantification of technetium-99 in wastewater by means of automated on-line extraction chromatography – anion-exchange chromatography – inductively coupled plasma-mass spectrometry »

J. Anal. At. Spectrom., 2024, 39, 2774-2782
Maximilian Horstmann, C. Derrick Quarles Jr, Steffen Happel, Michael Sperling, Andreas Faust, David Clases* and Uwe Karst

[1] S. Wagner, J. Santner, J. Irrgeher et al. Selective Diffusive Gradients in Thin Films (DGT) for the Simultaneous Assessment of Labile Sr and Pb Concentrations and Isotope Ratios in Soils. Anal Chem, 2022, 94, 16, 6338-6346, <https://doi.org/10.1021/acs.analchem.2c00546>

[2] S. Wagner, J. Santner, M. Puschenreiter, et al. Prediction of Radiogenic Sr and Pb Isotope Signatures in Plants using Diffusive Gradients in Thin Films. Analytical and Bioanalytical Chemistry, 2026, 418, 1541–1555, <https://doi.org/10.1007/s00216-026-06315-6>

[3] M. Horstmann, C. Derrick Quarles Jr, S. Happel et al. : Quantification of technetium-99 in wastewater by means of automated on-line extraction chromatography – anion-exchange chromatography – inductively coupled plasma-mass spectrometry. J. Anal. At. Spectrom., 2024, 39, 2774-2782, <https://doi.org/10.1039/D4JA00270A>





Development of alternative supports: impregnated discs already in use

Transposition of different resins extractive system onto discs to simplify preconcentration/purification/detection-counting

- Filtration or passive sampling on 25mm or 47mm
- Already in-use/published: TK100, TK201 and TK

Typical applications of TK-GA Discs [4]:

- pH 1 – 2 (for some actinides up to 3M HNO₃)
 - Disc preconditioning: 20% EtOH:H₂O, H₂O, then loading medium
 - Sample Loading, flow rate: 1 - 10mL/min (Lower flow rates favors better resolution)
 - 25 mm Discs => up to 100mL samples
 - 47 mm Discs => up to 1L (higher volumes possible depending on extractant system)
 - Rinsing with H₂O and 20% EtOH: H₂O
 - Drying and glueing onto steel disc => alpha spec
- Optional: subsequent Alpha/Beta by LSC
=> 100mL sample processed in 1,5h prior to alpha spec



[4] C. Bailly et al. Selective impregnated Membrane Filters for Direct Alpha Spectrometric Analysis of Actinides; Applied Radiation and Isot








Development of alternative supports:


Impregnated TLC (iTLC) supports for Fast Radionuclidic purity QC Checks

DGA Sheets by Kozempel et al. (CTUP) [5][6]

QC of radionuclides and generator eluents (e.g.: Ra-223, Ac-225/Bi-213, Pb-212,...)



I. $^{223}\text{Ra}[\text{RaCl}_2]$ migration ($R_f=0$)
 II. $^{223}\text{Ra}[\text{Ra}(\text{macropa})]$ ($R_f \approx 1$)
 III. ^{223}Ra tentatively chelated with DOTA showing unsuccessful radiolabeling (85% remained unchelated)
 IV. ^{223}Ra partially chelated with EDTA (65% remained unchelated). [5]



[5] D. S. Abou, N. A. Thiele, et al.: Towards the stable chelation of radium for biomedical applications with an 18-membered macrocyclic ligand. Chem. Sci., 2021, 12, 3733-3742, <https://doi.org/10.1039/D0SC06867E>

[6] J. Svedjehed et al.: "New extractant impregnated iTLC-SG paper facilitates improved TLC analysis for Cu radiolabelled peptides", poster presented at TERACHEM 2022, 14 – 17 September 2022, Bressanone (Italy). Poster available on <https://www.triskem-international.com/>



Development of alternative supports:

Impregnated TLC supports for Fast Radionuclidic purity QC Checks

CU Sheets by Svedjehed et al. [6]:

QC of Cu radiolabeled peptides (labeled vs free Cu): Both iTLC paper (SG and SG-CU) developed in less than 10 min. while Whatman paper took 25 – 30 min.

On-going studies on iTLC sheets : Zr, TK201, ... // Possible use of 2D TLC for radionuclide screening




Figure 1. TLC scans of [^{61}Cu]Cu-NOTA-octreotide spotted on: top, iTLC-SG; middle, Whatman paper; bottom extractant-impregnated iTLC-SG. Example 1 notes elevated levels of unlabelled ^{61}Cu , while example 2 notes comparable levels of labelled to unlabelled ^{61}Cu .

[6] J. Svedjehed et al.: "New extractant impregnated iTLC-SG paper facilitates improved TLC analysis for Cu radiolabelled peptides", poster presented at TERACHEM 2022, 14 – 17 September 2022, Bressanone (Italy). Poster available on <https://www.triskem-international.com/>



Extension to New Supports: impregnated filtering membranes

For impregnated membranes under developments (e.g. TK227):

➤ **See Inés Llopart presentation on applications**

On-going studies on extractive membranes :

- TK221 and TK200 for Actinides rapid screening
- new application using TK200 for Po in Ac-225
- TK101 for Ra and TK102 for Sr
- ...



Summary

- *Context of TK227 Resin development*
- *Characterisation of TK227 Resin*
- *Applications*
- *Development of alternative supports*
- **Conclusions**



Conclusions

The development of TK227 Resin has shown

- ✓ Y loading at 0,1-0,2M HNO_3 is possible
- ✓ Y is well separated from other elements
- ✓ Y needs to be carried
- It needs further testing (e.g. @ 0,5M-1M HNO_3) to avoid Y breakthrough in 1L up to 50L
- On-going methods developments for preconcentration/purification on resins and alternative supports

Thank you for your attention!

support@triskem.fr

